



plantmedia
a division of bioWORLD



1.888.bio.PLUS
www.bio-world.com



Ordering Information

Telephone Ordering

Call Toll-Free: 1.888.bio.PLUS (246-7587)
8:30 AM to 6:00 PM, Eastern-Standard Time
International: 1.614.792.8680

Order Online

www.bio-world.com **Fax Orders 1.614.792.8685**

Email Address

customerservice@bio-world.com

Mailing Address

4150 Tuller Rd. Suite 228, Dublin, OH 43017, USA

Technical Info

tech@bio-world.com

Shipping Charges

Oversize express service available for additional charge.

Delivery

In stock orders received Monday thru Friday before 2:00 PM are shipped the same day via Fedex or UPS ground. Next business day delivery within the carrier's overnight service area is available.

Discounts

Discount opportunities include the following: Receive a 10% discount when purchasing 10 units of any single item. Bulk (OEM) discounts and new lab discounts are also available, please inquire. A minimum order amount is \$50.00.

Terms

Net 30 days. Payments due in U.S. dollars. All international orders must be prepaid before shipment. For current prices, please visit www.bio-world.com.

Conditions

All products are sold for research or laboratory use only and not for humans or for any pharmaceutical purposes. Antibiotics are shipped only to institution/company laboratory addresses.

Research Use Only

Since these products are intended for research purposes by qualified persons, the Environmental Protection Agency does not require bioWorld to supply a premanufacturing notice.

Notice to Buyer/User

Nothing in this catalog is to be construed as an endorsement of any practice or product in violation of patents or regulations. It is the user's responsibility to determine for themselves the suitability of any material and/or procedure for a specific purpose and to adopt such safety precautions as may be necessary.

Warranty

bioWorld warrants that the product will meet the specifications stated on the technical data sheet. We agree to replace the product free of charge if the product does not conform to the specifications. This remedy is in lieu of all other remedies or claims for damages, consequential or otherwise, which the buyer may have against bioWorld.

Plant Tissue Culture

• transformation • regeneration • propagation

Products Overview

Plant cells are totipotent, ie., they are able to give rise to complete organism if provided with the necessary nutrients. Media used in plant tissue culture generally consist of mineral salts, sugars, vitamins, amino acids, growth regulators, and a semi-solid support of agar or gelrite. The minerals are needed by plant cells as building blocks for the synthesis of organic molecules, or as catalysts for enzymatic reactions. The ions in the media play an important role in the transport of molecules across the membranes, in osmotic regulation and in maintaining the electrochemical potential across the cell membranes.

Mineral Elements

The minerals in plant tissue culture media are divided into macro- and micro-elements. Ca, P, S, Mg, N and K are considered to be macro-elements while Fe, Cu, Mn, Co, Mo, B, I, Ni, Cl and Al are considered micro-elements. This categorization of elements into macro- and micro- is based on the quantities of the elements necessary to support healthy cell growth.

The quantity of micro-elements (or trace elements) in the media is not a representation of the importance of these elements to the plant. For example, an iron deficiency can have catastrophic effects on the growth and development of a plant. The necessity of micro-elements such as cobalt and aluminum is not very clear but certainly their presence in trace amounts in the medium display some growth benefits.

It is necessary that the plant cells receive all nutrients in balance thus making cations and anions equally available to support plant cell growth. For example, of potassium nitrate, both (K) and (NO) are essential. Within the group of macro-elements the nitrogen containing salts, such as potassium-, ammonium- and calcium nitrate, are among the most plentiful. Ammonium ions in excess may inhibit growth. All macro and micro elements must be supplied in a pure form, free of other elements and heavy metal contaminants.

Vitamins

Vitamins are needed in the plant media in small concentrations and must be supplied in natural forms. These compounds are essential for many biochemical reactions.

All media for plant cell and tissue culture contains thiamine (vitamin B1). Linsmaier and Skoog, after a thorough revision of the vitamins present in the MS medium, concluded that Thiamine is essential to plant growth.

Inositol significantly stimulates the growth and development of plants but it does not appear to be essential for growth. The effect of vitamins on the development of the cell in vitro differs from species to species, and therefore the concentration of vitamins needs to be adjusted for each explant.

Plantmedia provides freshly prepared vitamins for efficient growth of plant tissue in culture.

Call 1.888.BIOPLUS to order.

on • regeneration • propagation • transformation • regeneration • propagation

SKU	Product Description	Pathogen
Neogen Phytodetection		
30640000	ACLSV - Reagent set -	Apple Chlorotic Leafspot Virus
30640003	ACLSV - IDENTIKIT -	Apple Chlorotic Leafspot Virus
30640004	ACLSV - TOTAL Kit -	Apple Chlorotic Leafspot Virus
30640007	ACLSV - Positive control -	Apple Chlorotic Leafspot Virus
30640008	ACLSV - Negative control -	Apple Chlorotic Leafspot Virus
30640009	ACMV - Reagent set -	African Cassava Mosaic Virus
30640012	ACMV - IDENTIKIT -	African Cassava Mosaic Virus
30640013	ACMV - TOTAL Kit -	African Cassava Mosaic Virus
30640016	ACMV - Positive control -	African Cassava Mosaic Virus
30640017	ACMV - Negative control -	African Cassava Mosaic Virus
30640018	AgMV - Reagent set -	Agropyron Mosaic Virus
30640021	AgMV - IDENTIKIT -	Agropyron Mosaic Virus
30640022	AgMV - TOTAL Kit -	Agropyron Mosaic Virus
30640025	AgMV - Positive control -	Agropyron Mosaic Virus
30640026	AgMV - Negative control -	Agropyron Mosaic Virus
30640027	AMV - Reagent set -	Alfalfa mosaic virus
30640030	AMV - IDENTIKIT -	Alfalfa mosaic virus
30640031	AMV - TOTAL Kit -	Alfalfa mosaic virus
30640034	AMV - Positive control -	Alfalfa mosaic virus
30640035	AMV - Negative control -	Alfalfa mosaic virus
30640036	APLV - Reagent set -	Andean potato latent virus
30640039	APLV - IDENTIKIT -	Andean potato latent virus
30640040	APLV - TOTAL Kit -	Andean potato latent virus
30640043	APLV - Positive control -	Andean potato latent virus
30640044	APLV - Negative control -	Andean potato latent virus
30640045	APMoV - Reagent set -	Andean potato mottle virus
30640048	APMoV - IDENTIKIT -	Andean potato mottle virus
30640049	APMoV - TOTAL Kit -	Andean potato mottle virus
30640052	APMoV - Positive control -	Andean potato mottle virus
30640053	APMoV - Negative control -	Andean potato mottle virus
30640054	ApMV - Reagent set -	Apple mosaic virus
30640057	ApMV - IDENTIKIT -	Apple mosaic virus
30640058	ApMV - TOTAL Kit -	Apple mosaic virus
30640061	ApMV - Positive control -	Apple mosaic virus
30640062	ApMV - Negative control -	Apple mosaic virus
30640063	ArMV - Reagent Set -	Arabis mosaic virus
30640066	ArMV - IDENTIKIT -	Arabis mosaic virus
30640067	ArMV - TOTAL Kit -	Arabis mosaic virus
30640070	ArMV - Positive control -	Arabis mosaic virus
30640071	ArMV - Negative control -	Arabis mosaic virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640072	AVBO - Reagent set -	Arracacha virus B-OCA strain
30640075	AVBO - IDENTIKIT -	Arracacha virus B-OCA strain
30640076	AVBO - TOTAL Kit -	Arracacha virus B-OCA strain
30640079	AVBO - Positive control -	Arracacha virus B-OCA strain
30640080	AVBO - Negative control -	Arracacha virus B-OCA strain
30640081	BaMMV - Reagent set -	Barley mild mosiac virus
30640084	BaMMV - IDENTIKIT -	Barley mild mosiac virus
30640085	BaMMV - TOTAL Kit -	Barley mild mosiac virus
30640088	BaMMV - Positive control -	Barley mild mosiac virus
30640089	BaMMV - Negative control -	Barley mild mosiac virus
30640090	BaMMV - SPOTCHECK LF -	Barley mild mosiac virus
30640092	BaYMV - Reagent set -	Barley yellow mosiac virus
30640095	BaYMV - IDENTIKIT -	Barley yellow mosiac virus
30640096	BaYMV - TOTAL Kit -	Barley yellow mosiac virus
30640099	BaYMV - Positive control -	Barley yellow mosiac virus
30640100	BaYMV - Negative control -	Barley yellow mosiac virus
30640101	BBWV - Reagent set -	Broad bean wilt virus
30640104	BBWV - IDENTIKIT -	Broad bean wilt virus
30640105	BBWV - TOTAL Kit -	Broad bean wilt virus
30640108	BBWV - Positive control -	Broad bean wilt virus
30640109	BBWV - Negative control -	Broad bean wilt virus
30640110	BCMV - Reagent Set -	Bean common mosaic necrosis virus
30640113	BCMV - IDENTIKIT -	Bean common mosaic necrosis virus
30640114	BCMV - TOTAL Kit -	Bean common mosaic necrosis virus
30640117	BCMV - Positive control -	Bean common mosaic necrosis virus
30640118	BCMV - Negative control -	Bean common mosaic necrosis virus
30640119	BCMV - Reagent set -	Bean common mosaic virus
30640122	BCMV - IDENTIKIT -	Bean common mosaic virus
30640123	BCMV - TOTAL Kit -	Bean common mosaic virus
30640126	BCMV - Positive control -	Bean common mosaic virus
30640127	BCMV - Negative control -	Bean common mosaic virus
30640128	BGMV - Reagent set -	Bean golden mosaic virus
30640131	BGMV - IDENTIKIT -	Bean golden mosaic virus
30640132	BGMV - TOTAL Kit -	Bean golden mosaic virus
30640135	BGMV - Positive control -	Bean golden mosaic virus
30640136	BGMV - Negative control -	Bean golden mosaic virus

SKU	Product Description	Pathogen
30640137	BLRV - Reagent set -	Bean leafroll virus
30640140	BLRV - IDENTIKIT -	Bean leafroll virus
30640141	BLRV - TOTAL Kit -	Bean leafroll virus
30640144	BLRV - Positive control -	Bean leafroll virus
30640145	BLRV - Negative control -	Bean leafroll virus
30640146	BIShV - Reagent set -	Blueberry necrotic shock virus
30640149	BIShV - IDENTIKIT -	Blueberry necrotic shock virus
30640150	BIShV - TOTAL Kit -	Blueberry necrotic shock virus
30640153	BIShV - Positive control -	Blueberry necrotic shock virus
30640154	BIShV - Negative control -	Blueberry necrotic shock virus
30640155	BLSHV - Reagent set -	Blueberry shock virus
30640158	BLSHV - IDENTIKIT -	Blueberry shock virus
30640159	BLSHV - TOTAL Kit -	Blueberry shock virus
30640162	BLSHV - Positive control -	Blueberry shock virus
30640163	BLSHV - Negative control	Blueberry shock virus
30640164	BMYV - Reagent set-	Beet mild yellowing virus (syn. Beet western yellows - all strains)
30640167	BMYV - IDENTIKIT -	Beet mild yellowing virus (syn. Beet western yellows - all strains)
30640168	BMYV - TOTAL Kit -	Beet mild yellowing virus (syn. Beet western yellows - all strains)
30640171	BMYV - Positive control -	Beet mild yellowing virus (syn. Beet western yellows - all strains)
30640173	BNYVV - Reagent set -	Beet necrotic yellow vein virus
30640176	BNYVV - IDENTIKIT -	Beet necrotic yellow vein virus
30640177	BNYVV - TOTAL Kit -	Beet necrotic yellow vein virus
30640180	BNYVV - Positive control -	Beet necrotic yellow vein virus
30640181	BNYVV - Negative control -	Beet necrotic yellow vein virus
30640182	BNYVV - SPOTCHECK LF -	Beet necrotic yellow vein virus
30640184	BPMV - Reagent set -	Bean pod mottle virus
30640187	BPMV - IDENTIKIT -	Bean pod mottle virus
30640188	BPMV - TOTAL Kit -	Bean pod mottle virus
30640191	BPMV - Positive control -	Bean pod mottle virus
30640192	BPMV - Negative control -	Bean pod mottle virus
30640193	BSBV Reagent set -	Beet soil-borne virus
30640196	BSBV - IDENTIKIT -	Beet soil-borne virus
30640197	BSBV - TOTAL Kit -	Beet soil-borne virus
30640200	BSBV - Positive control -	Beet soil-borne virus
30640201	BSBV - Negative control -	Beet soil-borne virus
30640202	BScV - Reagent set -	Blueberry scorch virus
30640205	BScV - IDENTIKIT -	Blueberry scorch virus
30640206	BScV - TOTAL Kit -	Blueberry scorch virus
30640209	BScV - Positive control -	Blueberry scorch virus
30640210	BScV - Negative control -	Blueberry scorch virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640211	BSMV - Reagent set -	Barley stripe mosaic virus
30640214	BSMV - IDENTIKIT -	Barley stripe mosaic virus
30640215	BSMV - TOTAL Kit -	Barley stripe mosaic virus
30640218	BSMV - Positive control -	Barley stripe mosaic virus
30640219	BSMV - Negative control -	Barley stripe mosaic virus
30640220	BSSV - Reagent set -	Blueberry shoestring virus
30640223	BSSV - IDENTIKIT -	Blueberry shoestring virus
30640224	BSSV - TOTAL Kit -	Blueberry shoestring virus
30640227	BSSV - Positive control -	Blueberry shoestring virus
30640228	BSSV - Negative control -	Blueberry shoestring virus
30640229	BSV - Reagent set -	Banana streak virus
30640232	BSV - IDENTIKIT -	Banana streak virus
30640233	BSV - TOTAL Kit -	Banana streak virus
30640236	BSV - positive control -	Banana streak virus
30640237	BSV - Negative control -	Banana streak virus
30640238	BtMV - Reagent set -	Beet mosaic virus
30640241	BtMV - IDENTIKIT -	Beet mosaic virus
30640242	BtMV - TOTAL Kit -	Beet mosaic virus
30640245	BtMV - Positive control -	Beet mosaic virus
30640246	BtMV - Negative control -	Beet mosaic virus
30640247	BYDV - Screen Reagent set	Barley yellow dwarf virus screen
30640250	BYDV - Screen IDENTIKIT -	Barley yellow dwarf virus screen
30640251	BYDV - Screen TOTAL Kit -	Barley yellow dwarf virus screen
30640254	BYDV - Screen Positive control -	Barley yellow dwarf virus screen
30640255	BYDV - Screen Negative control -	Barley yellow dwarf virus screen
30640256	BYDV-MAV - Reagent set	Barley yellow dwarf virus - MAV
30640259	BYDV-MAV - IDENTIKIT -	Barley yellow dwarf virus - MAV
30640260	BYDV-MAV - TOTAL Kit -	Barley yellow dwarf virus - MAV
30640263	BYDV-MAV - Positive control -	Barley yellow dwarf virus - MAV
30640264	BYDV-MAV - Negative control -	Barley yellow dwarf virus - MAV
30640265	BYDV-PAV - Reagent set -	Barley yellow dwarf virus - PAV
30640268	BYDV-PAV - IDENTIKIT -	Barley yellow dwarf virus - PAV
30640269	BYDV-PAV - TOTAL Kit -	Barley yellow dwarf virus - PAV
30640272	BYDV-PAV - Positive control -	Barley yellow dwarf virus - PAV
30640273	BYDV-PAV - Negative control -	Barley yellow dwarf virus - PAV
30640274	BYDV-RPV - Reagent set -	Barley yellow dwarf virus - RPV
30640277	BYDV-RPV - IDENTIKIT -	Barley yellow dwarf virus - RPV
30640278	BYDV-RPV - TOTAL Kit -	Barley yellow dwarf virus - RPV
30640281	BYDV-RPV - Positive control -	Barley yellow dwarf virus - RPV
30640282	BYDV-RPV - Negative control -	Barley yellow dwarf virus - RPV

Call 1.888.BIOPRO to order.

SKU	Product Description	Pathogen
30640283	BYMV - Reagent set -	Bean yellow mosaic virus
30640286	BYMV - Identikit -	Bean yellow mosaic virus
30640287	BYMV - TOTAL Kit -	Bean yellow mosaic virus
30640290	BYMV - Positive control -	Bean yellow mosaic virus
30640291	BYMV - Negative control -	Bean yellow mosaic virus
30640292	BYV - Reagent Set -	Beet yellow virus
30640295	BYV - IDENTIKIT -	Beet yellow virus
30640296	BYV - TOTAL Kit -	Beet yellow virus
30640299	BYV - Positive control -	Beet yellow virus
30640300	BYV - Negative control -	Beet yellow virus
30640301	CaMV - Reagent set-	Cauliflower mosaic virus
30640304	CaMV - IDENTIKIT -	Cauliflower mosaic virus
30640305	CaMV - TOTAL Kit -	Cauliflower mosaic virus
30640308	CaMV - Positive control -	Cauliflower mosaic virus
30640309	CaMV - Negative control -	Cauliflower mosaic virus
30640310	CarMV - Reagent set -	Carnation mottle virus
30640313	CarMV - IDENTIKIT -	Carnation mottle virus
30640314	CarMV - TOTAL Kit -	Carnation mottle virus
30640317	CarMV - Positive control -	Carnation mottle virus
30640318	CarMV - Negative control -	Carnation mottle virus
30640319	CarVY - Reagent set -	Carrot virus Y
30640322	CarVY - IDENTIKIT -	Carrot virus Y
30640323	CarVY - TOTAL Kit -	Carrot virus Y
30640326	CarVY - Positive control -	Carrot virus Y
30640327	CarVY - Negative control -	Carrot virus Y
30640328	CERV - Reagent Set -	Carnation etched ring virus
30640331	CERV - INDENTIKIT -	Carnation etched ring virus
30640332	CERV - TOTAL Kit -	Carnation etched ring virus
30640335	CERV - Positive control -	Carnation etched ring virus
30640336	CERV - Negative control -	Carnation etched ring virus
30640337	CGMMV - Reagent set -	Cucumber green mottle mosaic virus
30640340	CGMMV - IDENTIKIT -	Cucumber green mottle mosaic virus
30640341	CGMMV - TOTAL KIT -	Cucumber green mottle mosaic virus
30640344	CGMMV - Positive control -	Cucumber green mottle mosaic virus
30640345	CGMMV - Negative control -	Cucumber green mottle mosaic virus
30640346	CGMMV - SPOTCHECK LF -	Cucumber green mottle mosaic virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640348	ChiVMV - Reagent set -	Chilli veinal mottle virus
30640351	ChiVMV - IDENTIKIT -	Chilli veinal mottle virus
30640352	ChiVMV - TOTAL Kit -	Chilli veinal mottle virus
30640355	ChiVMV - Positive control -	Chilli veinal mottle virus
30640356	ChiVMV - Negative control -	Chilli veinal mottle virus
30640357	CLRV - Reagent Set -	Cherry leaf roll virus
30640360	CLRV - IDENTIKIT -	Cherry leaf roll virus
30640361	CLRV - TOTAL Kit -	Cherry leaf roll virus
30640364	CLRV - Positive control -	Cherry leaf roll virus
30640365	CLRV - Negative control -	Cherry leaf roll virus
30640366	CLV - Reagent set -	Carnation latent virus
30640369	CLV - IDENTIKIT -	Carnation latent virus
30640370	CLV - TOTAL Kit -	Carnation latent virus
30640373	CLV - Positive control -	Carnation latent virus
30640374	CLV - Negative control -	Carnation latent virus
30640375	CMV - Reagent set -	Cucumber mosaic virus
30640378	CMV - IDENTIKIT -	Cucumber mosaic virus
30640379	CMV - TOTAL Kit -	Cucumber mosaic virus
30640382	CMV - Positive control -	Cucumber mosaic virus
30640383	CMV - Negative control -	Cucumber mosaic virus
30640384	CMV - SPOTCHECK LF -	Cucumber mosaic virus
30640386	CNFV - Reagent set -	Carnation necrotic fleck virus
30640389	CNFV - IDENTIKIT -	Carnation necrotic fleck virus
30640390	CNFV - TOTAL Kit -	Carnation necrotic fleck virus
30640393	CNFV - Positive control -	Carnation necrotic fleck virus
30640394	CNFV - Negative control -	Carnation necrotic fleck virus
30640395	CpABMV - Reagent Set -	Cowpea aphid-borne mosaic virus
30640398	CpABMV - IDENTIKIT -	Cowpea aphid-borne mosaic virus
30640399	CpABMV - TOTAL Kit -	Cowpea aphid-borne mosaic virus
30640402	CpABMV - Positive control -	Cowpea aphid-borne mosaic virus
30640403	CpABMV - Negative control -	Cowpea aphid-borne mosaic virus
30640404	CpCMoV - Reagent set -	Cowpea chlorotic mottle virus
30640407	CpCMoV - IDENTIKIT -	Cowpea chlorotic mottle virus
30640408	CpCMoV - TOTAL Kit -	Cowpea chlorotic mottle virus
30640411	CpCMoV - Positive control -	Cowpea chlorotic mottle virus
30640412	CpCMoV - Negative control -	Cowpea chlorotic mottle virus

SKU	Product Description	Pathogen
30640413	CpMMov - Reagent Set -	Cowpea mild mottle virus
30640416	CpMMoV - IDENTIKIT -	Cowpea mild mottle virus
30640417	CpMMoV - TOTAL Kit -	Cowpea mild mottle virus
30640420	CPMMoV - Positive control -	Cowpea mild mottle virus
30640421	CpMMoV - Negative control -	Cowpea mild mottle virus
30640422	CpMV - Reagent set -	Cowpea mosaic virus
30640425	CpMV - IDENTIKIT -	Cowpea mosaic virus
30640426	CpMV - TOTAL Kit -	Cowpea mosaic virus
30640429	CpMV - Positive control -	Cowpea mosaic virus
30640430	CpMV - Negative control -	Cowpea mosaic virus
30640431	CpSMV - Reagent Set -	Cowpea severe mosaic virus
30640434	CPSMV - IDENTIKIT -	Cowpea severe mosaic virus
30640435	CpSMV - TOTAL Kit -	Cowpea severe mosaic virus
30640438	CpSMV - Positive control -	Cowpea severe mosaic virus
30640439	CpSMV - Negative control -	Cowpea severe mosaic virus
30640440	CRLV - Reagent set -	Cherry rasp leaf virus
30640443	CRLV - IDENTIKIT -	Cherry rasp leaf virus
30640444	CRLV - TOTAL Kit -	Cherry rasp leaf virus
30640447	CRLV - Positive control -	Cherry rasp leaf virus
30640448	CRLV - Negative control -	Cherry rasp leaf virus
30640449	CRSV - Reagent Set -	Carnation ringspot virus
30640452	CRSV - IDENTIKIT -	Carnation ringspot virus
30640453	CRSV - TOTAL Kit -	Carnation ringspot virus
30640456	CRSV - Positive control -	Carnation ringspot virus
30640457	CRSV - Negative control -	Carnation ringspot virus
30640458	CSNV - Reagent set -	Chrysanthemum stem necrosis virus
30640461	CSNV - IDENTIKIT -	Chrysanthemum stem necrosis virus
30640462	CSNV - TOTAL Kit -	Chrysanthemum stem necrosis virus
30640465	CSNV - Positive control -	Chrysanthemum stem necrosis virus
30640466	CSNV - Negative control -	Chrysanthemum stem necrosis virus
30640467	CTV - Reagent set -	Citrus tristeza virus
30640470	CTV - IDENTIKIT -	Citrus tristeza virus
30640471	CTV - TOTAL Kit -	Citrus tristeza virus
30640474	CTV - Positive control -	Citrus tristeza virus
30640475	CTV - Negative control -	Citrus tristeza virus
30640476	CVB - Reagent set -	Chrysanthemum virus B
30640479	CVB - IDENTIKIT -	Chrysanthemum virus B
30640480	CVB - TOTAL Kit -	Chrysanthemum virus B
30640483	CVB - Positive control -	Chrysanthemum virus B
30640484	CVB - Negative control -	Chrysanthemum virus B

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640485	CVMV - Reagent set -	Carnation vein mottle virus
30640488	CVMV - IDENTIKIT -	Carnation vein mottle virus
30640489	CVMV - TOTAL Kit -	Carnation vein mottle virus
30640492	CVMV - Positive control -	Carnation vein mottle virus
30640493	CVMV - Negative control -	Carnation vein mottle virus
30640494	CymMV - Reagent set -	Cymbidium mosaic virus
30640497	CymMV - IDENTIKIT -	Cymbidium mosaic virus
30640498	CymMV - TOTAL Kit -	Cymbidium mosaic virus
30640501	CymMV - Positive control -	Cymbidium mosaic virus
30640502	CymMV - Negative control -	Cymbidium mosaic virus
30640503	CymMV - SPOTCHECK LF -	Cymbidium mosaic virus
30640505	CYVV - Reagent set -	Clover yellow vein virus
30640508	CYVV - IDENTIKIT -	Clover yellow vein virus
30640509	CYVV - TOTAL Kit -	Clover yellow vein virus
30640512	CYVV - Positive control -	Clover yellow vein virus
30640513	CYVV - Negative control -	Clover yellow vein virus
30640514	EPMV - Reagent set -	Egg plant mosaic virus
30640517	EPMV - IDENTIKIT -	Egg plant mosaic virus
30640518	EPMV - TOTAL Kit -	Egg plant mosaic virus
30640521	EPMV - Positive control -	Egg plant mosaic virus
30640522	EPMV - Negative control -	Egg plant mosaic virus
30640532	GarCLV - Reagent set -	Garlic common latent virus
30640535	GarCLV - IDENTIKIT -	Garlic common latent virus
30640536	GarCLV - TOTAL Kit -	Garlic common latent virus
30640539	GarCLV - Positive control -	Garlic common latent virus
30640540	GarCLV - Negative control -	Garlic common latent virus
30640541	Gemini Virus Screen - Reagent set -	Gemini virus screen
30640544	Gemini Virus Screen - IDENTIKIT -	Gemini virus screen
30640545	Gemini Virus Screen - TOTAL Kit -	Gemini virus screen
30640548	Gemini Virus Screen - Positive control -	Gemini virus screen
30640549	Gemini Virus Screen - Negative control -	Gemini virus screen
30640550	GFKV - Reagent set -	Grapevine fleck virus
30640553	GFKV - IDENTIKIT -	Grapevine fleck virus
30640554	GFKV - TOTAL Kit -	Grapevine fleck virus
30640557	GFKV - Positive control -	Grapevine fleck virus
30640558	GFKV - Negative control -	Grapevine fleck virus

SKU	Product Description	Pathogen
30640559	GFLV - Reagent set -	Grapevine fanleaf virus
30640562	GFLV - IDENTIKIT -	Grapevine fanleaf virus
30640563	GFLV - TOTAL Kit -	Grapevine fanleaf virus
30640566	GFLV - Positive control -	Grapevine fanleaf virus
30640567	GFLV - Negative control -	Grapevine fanleaf virus
30640568	GLRaV-1 - Reagent set -	Grapevine leafroll virus 1
30640571	GLRaV-1 - IDENTIKIT -	Grapevine leafroll virus 1
30640572	GLRaV-1 - TOTAL Kit -	Grapevine leafroll virus 1
30640575	GLRaV-1 - Positive control -	Grapevine leafroll virus 1
30640576	GLRaV-1 - Negative control -	Grapevine leafroll virus 1
30640577	GLRaV-1,3 - Reagent set -	Grapevine leafroll virus 1, 3
30640580	GLRaV-1,3 - IDENTIKIT -	Grapevine leafroll virus 1, 3
30640581	GLRaV-1,3 - TOTAL Kit -	Grapevine leafroll virus 1, 3
30640584	GLRaV-1,3 - Positive control -	Grapevine leafroll virus 1, 3
30640585	GLRaV-1,3 - Negative control -	Grapevine leafroll virus 1, 3
30640586	GLRaV-2 - Reagent set -	Grapevine leafroll virus 2
30640589	GLRaV-2 - IDENTIKIT -	Grapevine leafroll virus 2
30640590	GLRaV-2 - TOTAL Kit -	Grapevine leafroll virus 2
30640593	GLRaV-2 - Positive control -	Grapevine leafroll virus 2
30640594	GLRaV-2 - Negative control -	Grapevine leafroll virus 2
30640595	GLRaV-3 - Reagent set -	Grapevine leafroll virus 3
30640598	GLRaV-3 - IDENTIKIT -	Grapevine leafroll virus 3
30640599	GLRaV-3 - TOTAL Kit -	Grapevine leafroll virus 3
30640602	GLRaV-3 - Positive control -	Grapevine leafroll virus 3
30640603	GLRaV-3 - Negative control -	Grapevine leafroll virus 3
30640604	GLRaV-5 - Reagent set -	Grapevine leafroll virus 5
30640607	GLRaV-5 - IDENTIKIT -	Grapevine leafroll virus 5
30640608	GLRaV-5 - TOTAL Kit -	Grapevine leafroll virus 5
30640611	GLRaV-5 - Positive control -	Grapevine leafroll virus 5
30640612	GLRaV-5 - Negative control -	Grapevine leafroll virus 5
30640613	GLRaV-7 - Reagent set -	Grapevine leafroll virus 7
30640616	GLRaV-7 - IDENTIKIT -	Grapevine leafroll virus 7
30640617	GLRaV-7 - TOTAL Kit -	Grapevine leafroll virus 7
30640620	GLRaV-7 - Positive control -	Grapevine leafroll virus 7
30640621	GLRaV-7 - Negative control -	Grapevine leafroll virus 7
30640622	GVA - Reagent set -	Grapevine virus A
30640625	GVA - IDENTIKIT -	Grapevine virus A
30640626	GVA - TOTAL Kit -	Grapevine virus A
30640629	GVA - Positive control -	Grapevine virus A
30640630	GVA - Negative control -	Grapevine virus A

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640631	GVB - Reagent set -	Grapevine virus B
30640634	GVB - IDENTIKIT -	Grapevine virus B
30640635	GVB - TOTAL Kit -	Grapevine virus B
30640638	GVB - Positive control -	Grapevine virus B
30640639	GVB - Negative control -	Grapevine virus B
30640640	HVX - Reagent set -	Hosta virus X
30640643	HVX - IDENTIKIT -	Hosta virus X
30640644	HVX - TOTAL Kit -	Hosta virus X
30640647	HVX - Positive control -	Hosta virus X
30640648	HVX - Negative control -	Hosta virus X
30640649	INSV - Reagent set -	Impatiens necrotic spot virus
30640652	INSV - IDENTIKIT -	Impatiens necrotic spot virus
30640653	INSV - TOTAL Kit -	Impatiens necrotic spot virus
30640656	INSV - Positive control -	Impatiens necrotic spot virus
30640657	INSV - Negative control -	Impatiens necrotic spot virus
30640658	INSV - SPOTCHECK LF -	Impatiens necrotic spot virus
30640660	IYSV - Reagent set -	Iris yellow spot virus
30640663	IYSV - IDENTIKIT -	Iris yellow spot virus
30640664	IYSV - TOTAL Kit -	Iris yellow spot virus
30640667	IYSV - Positive control -	Iris yellow spot virus
30640668	IYSV - Negative control -	Iris yellow spot virus
30640669	KGMMV - Reagent set -	Kyuri green mottle mosaic virus
30640672	KGMMV - IDENTIKIT -	Kyuri green mottle mosaic virus
30640673	KGMMV - TOTAL Kit -	Kyuri green mottle mosaic virus
30640676	KGMMV - Positive control -	Kyuri green mottle mosaic virus
30640677	KGMMV - Negative control -	Kyuri green mottle mosaic virus
30640678	LMV - Reagent set -	Lettuce mosaic virus
30640681	LMV - IDENTIKIT -	Lettuce mosaic virus
30640682	LMV - TOTAL Kit -	Lettuce mosaic virus
30640685	LMV - Positive control -	Lettuce mosaic virus
30640686	LMV - Negative control -	Lettuce mosaic virus
30640687	LSV - Reagent set -	Lily symptomless virus
30640690	LSV - IDENTIKIT -	Lily symptomless virus
30640691	LSV - TOTAL Kit -	Lily symptomless virus
30640694	LSV - Positive control -	Lily symptomless virus
30640695	LSV - Negative control -	Lily symptomless virus

SKU	Product Description	Pathogen
30640696	LVX - Reagent set -	Lily virus X
30640699	LVX - IDENTIKIT -	Lily virus X
30640700	LVX - TOTAL Kit -	Lily virus X
30640703	LVX - Positive control -	Lily virus X
30640704	LVX - Negative control -	Lily virus X
30640705	LYSV - Reagent set -	Leek yellow stripe virus
30640708	LYSV - IDENTIKIT -	Leek yellow stripe virus
30640709	LYSV - TOTAL Kit -	Leek yellow stripe virus
30640712	LYSV - Positive control -	Leek yellow stripe virus
30640713	LYSV - Negative control -	Leek yellow stripe virus
30640714	MCDV - Reagent set -	Maize chlorotic dwarf virus
30640717	MCDV - IDENTIKIT -	Maize chlorotic dwarf virus
30640718	MCDV - TOTAL Kit -	Maize chlorotic dwarf virus
30640721	MCDV - Positive control -	Maize chlorotic dwarf virus
30640722	MCDV - Negative control -	Maize chlorotic dwarf virus
30640723	MCMV - Reagent set -	Maize chlorotic mottle virus
30640726	MCMV - IDENTIKIT -	Maize chlorotic mottle virus
30640727	MCMV - TOTAL Kit -	Maize chlorotic mottle virus
30640730	MCMV - Positive control -	Maize chlorotic mottle virus
30640731	MCMV - Negative control -	Maize chlorotic mottle virus
30640732	MDMV - Reagent set -	Maize dwarf mosaic virus
30640735	MDMV - IDENTIKIT -	Maize dwarf mosaic virus
30640736	MDMV - TOTAL Kit -	Maize dwarf mosaic virus
30640739	MDMV - Positive control -	Maize dwarf mosaic virus
30640740	MDMV - Negative control -	Maize dwarf mosaic virus
30640741	MiLBVV - Reagent set -	Mirafiori lettuce big vein virus
30640744	MiLBVV - IDENTIKIT -	Mirafiori lettuce big vein virus
30640745	MiLBVV - TOTAL Kit -	Mirafiori lettuce big vein virus
30640748	MiLBVV - Positive control -	Mirafiori lettuce big vein virus
30640749	MiLBVV - Negative control -	Mirafiori lettuce big vein virus
30640750	MNSV - Reagent set -	Melon necrotic spot virus
30640753	MNSV - IDENTIKIT -	Melon necrotic spot virus
30640754	MNSV - TOTAL Kit -	Melon necrotic spot virus
30640757	MNSV - Positive control -	Melon necrotic spot virus
30640758	MNSV - Negative control -	Melon necrotic spot virus
30640759	MSV - Reagent set -	Maize streak virus
30640762	MSV - IDENTIKIT -	Maize streak virus
30640763	MSV - TOTAL Kit -	Maize streak virus
30640766	MSV - Positive control -	Maize streak virus
30640767	MSV - Negative control -	Maize streak virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640768	MV1 - Reagent set -	Mushroom virus 1
30640771	MV1 - IDENTIKIT -	Mushroom virus 1
30640772	MV1 - TOTAL Kit -	Mushroom virus 1
30640775	MV1 - Positive control -	Mushroom virus 1
30640776	MV1 - Negative control -	Mushroom virus 1
30640777	MWMV - Reagent set -	Moroccan watermelon mosaic virus
30640780	MWMV - IDENTIKIT -	Moroccan watermelon mosaic virus
30640781	MWMV - TOTAL Kit -	Moroccan watermelon mosaic virus
30640784	MWMV - Positive control -	Moroccan watermelon mosaic virus
30640785	MWMV - Negative control -	Moroccan watermelon mosaic virus
30640786	NLSYV - Reagent set -	Narcissus late season yellows virus
30640789	NLSYV - IDENTIKIT -	Narcissus late season yellows virus
30640790	NLSYV - TOTAL Kit -	Narcissus late season yellows virus
30640793	NLSYV - Positive control -	Narcissus late season yellows virus
30640794	NLSYV - Negative control -	Narcissus late season yellows virus
30640795	NLV - Reagent set -	Narcissus latent virus
30640798	NLV - IDENTIKIT -	Narcissus latent virus
30640799	NLV - TOTAL Kit -	Narcissus latent virus
30640802	NLV - Positive control -	Narcissus latent virus
30640803	NLV - Negative control -	Narcissus latent virus
30640804	NMV - Reagent set -	Narcissus mosaic virus
30640807	NMV - IDENTIKIT -	Narcissus mosaic virus
30640808	NMV - TOTAL Kit -	Narcissus mosaic virus
30640811	NMV - Positive control -	Narcissus mosaic virus
30640812	NMV - Negative control -	Narcissus mosaic virus
30640813	NQV - Reagent set -	Narcissus Q virus
30640816	NQV - IDENTIKIT -	Narcissus Q virus
30640817	NQV - TOTAL Kit -	Narcissus Q virus
30640820	NQV - Positive control -	Narcissus Q virus
30640821	NQV - Negative control -	Narcissus Q virus
30640822	NTNV - Reagent set -	Narcissus tip necrosis virus
30640825	NTNV - IDENTIKIT -	Narcissus tip necrosis virus
30640826	NTNV - TOTAL Kit -	Narcissus tip necrosis virus
30640829	NTNV - Positive control -	Narcissus tip necrosis virus
30640830	NTNV - Negative control -	Narcissus tip necrosis virus
30640831	NYSV - Reagent set -	Narcissus yellow stripe virus
30640834	NYSV - IDENTIKIT -	Narcissus yellow stripe virus
30640835	NYSV - TOTAL Kit -	Narcissus yellow stripe virus
30640838	NYSV - Positive control -	Narcissus yellow stripe virus
30640839	NYSV - Negative control -	Narcissus yellow stripe virus

Call 1.888.BIOPLUS to order.

SKU	Product Description	Pathogen
30640842	ORSV - Reagent set -	Odontoglossum ringspot virus
30640843	ORSV - IDENTIKIT -	Odontoglossum ringspot virus
30640844	ORSV - TOTAL Kit -	Odontoglossum ringspot virus
30640847	ORSV - Positive control -	Odontoglossum ringspot virus
30640848	ORSV - Negative control -	Odontoglossum ringspot virus
30640849	ORSV - SPOTCHECK LF -	Odontoglossum ringspot virus
30640851	OYDV - Reagent set -	Onion yellow dwarf virus
30640854	OYDV - IDENTIKIT -	Onion yellow dwarf virus
30640855	OYDV - TOTAL Kit -	Onion yellow dwarf virus
30640858	OYDV - Positive control -	Onion yellow dwarf virus
30640859	OYDV - Negative control -	Onion yellow dwarf virus
30640860	PAMV - Reagent set -	Potato acuba mosaic virus
30640863	PAMV - IDENTIKIT -	Potato acuba mosaic virus
30640864	PAMV - TOTAL Kit -	Potato acuba mosaic virus
30640867	PAMV - Positive control -	Potato acuba mosaic virus
30640868	PAMV - Negative control -	Potato acuba mosaic virus
30640869	PapMV - Reagent set -	Papaya mosaic virus
30640872	PapMV - IDENTIKIT -	Papaya mosaic virus
30640873	PapMV - TOTAL Kit -	Papaya mosaic virus
30640876	PapMV - Positive control -	Papaya mosaic virus
30640877	PapMV - Negative control -	Papaya mosaic virus
30640878	PBRSV - Reagent set -	Potato black ringspot virus
30640881	PBRSV - IDENTIKIT -	Potato black ringspot virus
30640882	PBRSV - TOTAL Kit -	Potato black ringspot virus
30640885	PBRSV - Positive control -	Potato black ringspot virus
30640886	PBRSV - Negative control -	Potato black ringspot virus
30640887	PDV - Reagent set -	Prune dwarf virus
30640890	PDV - IDENTIKIT -	Prune dwarf virus
30640891	PDV - TOTAL Kit -	Prune dwarf virus
30640894	PDV - Positive control -	Prune dwarf virus
30640895	PDV - Negative control -	Prune dwarf virus
30640896	PEMV - Reagent set -	Pea enation mosaic virus
30640899	PEMV - IDENTIKIT -	Pea enation mosaic virus
30640900	PEMV - TOTAL Kit -	Pea enation mosaic virus
30640903	PEMV - Positive control -	Pea enation mosaic virus
30640904	PEMV - Negative control -	Pea enation mosaic virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30640905	PepMoV - Reagent set -	Pepper mottle virus
30640908	PepMoV - IDENTIKIT -	Pepper mottle virus
30640909	PepMoV - TOTAL Kit -	Pepper mottle virus
30640912	PepMoV - Positive control -	Pepper mottle virus
30640913	PepMoV - Negative control -	Pepper mottle virus
30640914	PepMV - Reagent set -	Pepino mosaic virus
30640917	PepMV - IDENTIKIT -	Pepino mosaic virus
30640918	PepMV - TOTAL Kit -	Pepino mosaic virus
30640921	PepMV - Positive control -	Pepino mosaic virus
30640922	PepMV - Negative control -	Pepino mosaic virus
30640923	PepMV - SPOTCHECK LF -	Pepino mosaic virus
30640925	PFBV - Reagent set -	Pelargonium flower break virus
30640928	PFBV - IDENTIKIT -	Pelargonium flower break virus
30640929	PFBV - TOTAL Kit -	Pelargonium flower break virus
30640932	PFBV - Positive control -	Pelargonium flower break virus
30640933	PFBV - Negative control -	Pelargonium flower break virus
30640934	PLCV - Reagent set -	Pelargonium leaf curl virus
30640937	PLCV - IDENTIKIT -	Pelargonium leaf curl virus
30640938	PLCV - TOTAL Kit -	Pelargonium leaf curl virus
30640941	PLCV - Positive control -	Pelargonium leaf curl virus
30640942	PLCV - Negative control -	Pelargonium leaf curl virus
30640943	PLPV - Reagent set -	Pelargonium line pattern virus
30640946	PLPV - IDENTIKIT -	Pelargonium line pattern virus
30640947	PLPV - TOTAL Kit -	Pelargonium line pattern virus
30640950	PLPV - Positive control -	Pelargonium line pattern virus
30640951	PLPV - Negative control -	Pelargonium line pattern virus
30640952	PLRV - Reagent set -	Potato leafroll virus - DAS
30640956	PLRV - IDENTIKIT -	Potato leafroll virus - DAS
30640957	PLRV - TOTAL Kit -	Potato leafroll virus - DAS
30640957	PLRV - TOTAL Kit -	Potato leafroll virus - DAS
30640957	PLRV - TOTAL Kit -	Potato leafroll virus - DAS
30640960	PLRV - Positive control -	Potato leafroll virus - DAS
30640961	PLRV - Negative control -	Potato leafroll virus - DAS
30640962	PLRV - Reagent set -	Potato leafroll virus - TAS
30640966	PLRV - IDENTIKIT -	Potato leafroll virus - TAS
30640967	PLRV - TOTAL Kit -	Potato leafroll virus - TAS
30640970	PLRV - Positive control -	Potato leafroll virus - TAS
30640971	PLRV - Negative control -	Potato leafroll virus - TAS

SKU	Product Description	Pathogen
30640972	PLRV - Reagent set -	Potato leafroll virus
30640976	PLRV - IDENTIKIT -	Potato leafroll virus
30640977	PLRV - TOTAL Kit -	Potato leafroll virus
30640980	PLRV - Positive control -	Potato leafroll virus
30640981	PLRV - Negative control -	Potato leafroll virus
30640982	PMMoV - Reagent set -	Pepper mild mottle virus
30640985	PMMoV - IDENTIKIT -	Pepper mild mottle virus
30640986	PMMoV - TOTAL Kit -	Pepper mild mottle virus
30640989	PMMoV - Positive control -	Pepper mild mottle virus
30640990	PMMoV - Negative control -	Pepper mild mottle virus
30640991	PMTV - Reagent set -	Potato mop top virus - TAS
30640994	PMTV - IDENTIKIT -	Potato mop top virus - TAS
30640995	PMTV - TOTAL Kit -	Potato mop top virus - TAS
30640998	PMTV - Positive control -	Potato mop top virus - TAS
30640999	PMTV - Negative control -	Potato mop top virus - TAS
30641000	PMTV - Reagent set -	Potato mop top virus - DAS
30641003	PMTV - IDENTIKIT -	Potato mop top virus - DAS
30641004	PMTV - TOTAL Kit -	Potato mop top virus - DAS
30641007	PMTV - Positive control -	Potato mop top virus - DAS
30641008	PMTV - Negative control -	Potato mop top virus - DAS
30641009	PNRSV - Reagent set -	Prunus necrotic ringspot virus
30641012	PNRSV - IDENTIKIT -	Prunus necrotic ringspot virus
30641013	PNRSV - TOTAL Kit -	Prunus necrotic ringspot virus
30641016	PNRSV - Positive control -	Prunus necrotic ringspot virus
30641017	PNRSV - Negative control -	Prunus necrotic ringspot virus
30641018	Potyvirus - Reagent set -	Potyvirus Detector
30641021	Potyvirus - IDENTIKIT -	Potyvirus Detector
30641022	Potyvirus - TOTAL Kit -	Potyvirus Detector
30641025	Potyvirus - Positive control -	Potyvirus Detector
30641026	Potyvirus - Negative control -	Potyvirus Detector
30641027	PPV - Reagent set -	Plum pox virus
30641030	PPV - IDENTIKIT -	Plum pox virus
30641031	PPV - TOTAL Kit -	Plum pox virus
30641034	PPV - Positive control -	Plum pox virus
30641035	PPV - Negative control -	Plum pox virus
30641036	PPV - SPOTCHECK LF -	Plum pox virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641038	PRMV - Reagent set -	Peach rosette mosaic virus
30641041	PRMV - IDENTIKIT -	Peach rosette mosaic virus
30641042	PRMV - TOTAL Kit -	Peach rosette mosaic virus
30641045	PRMV - Positive control -	Peach rosette mosaic virus
30641046	PRMV - Negative control -	Peach rosette mosaic virus
30641047	PRSV - Reagent set -	Papaya ringspot virus
30641050	PRSV - IDENTIKIT -	Papaya ringspot virus
30641051	PRSV - TOTAL Kit -	Papaya ringspot virus
30641054	PRSV - Positive control -	Papaya ringspot virus
30641055	PRSV - Negative control -	Papaya ringspot virus
30641056	PSbMV - Reagent set -	Pea seed-borne mosaic virus
30641059	PSbMV - IDENTIKIT -	Pea seed-borne mosaic virus
30641060	PSbMV - TOTAL Kit -	Pea seed-borne mosaic virus
30641063	PSbMV - Positive control -	Pea seed-borne mosaic virus
30641064	PSbMV - Negative control -	Pea seed-borne mosaic virus
30641065	PSV - Reagent set -	Peanut stripe virus
30641068	PSV - IDENTIKIT -	Peanut stripe virus
30641069	PSV - TOTAL Kit -	Peanut stripe virus
30641072	PSV - Positive control -	Peanut stripe virus
30641073	PSV - Negative control -	Peanut stripe virus
30641074	PSV - Reagent set -	Peanut stunt virus
30641077	PSV - IDENTIKIT -	Peanut stunt virus
30641078	PSV - TOTAL Kit - 250	Peanut stunt virus
30641081	PSV - Positive control -	Peanut stunt virus
30641082	PSV - Negative control -	Peanut stunt virus
30641083	PVA - Reagent set -	Potato virus A
30641087	PVA - IDENTIKIT -	Potato virus A
30641088	PVA - TOTAL Kit -	Potato virus A
30641091	PVA - Positive control -	Potato virus A
30641092	PVA - Negative control -	Potato virus A
30641093	PVA - Reagent set -	Potato virus A
30641097	PVA - IDENTIKIT -	Potato virus A
30641098	PVA - TOTAL Kit -	Potato virus A
30641101	PVA - Positive control -	Potato virus A
30641102	PVA - Negative control -	Potato virus A
30641103	PVA - SPOTCHECK LF -	Potato virus A

SKU	Product Description	Pathogen
30641105	PVM - Reagent set -	Potato virus M
30641109	PVM - IDENTIKIT -	Potato virus M
30641110	PVM - TOTAL Kit -	Potato virus M
30641113	PVM - Positive control -	Potato virus M
30641114	PVM - Negative control -	Potato virus M
30641115	PVM - Reagent set -	Potato virus M
30641119	PVM - IDENTIKIT -	Potato virus M
30641120	PVM - TOTAL Kit -	Potato virus M
30641123	PVM - Positive control -	Potato virus M
30641124	PVM - Negative control -	Potato virus M
30641125	PVMoV - Reagent set -	Pepper veinal mottle virus
30641128	PVMoV - IDENTIKIT -	Pepper veinal mottle virus
30641129	PVMoV - TOTAL Kit -	Pepper veinal mottle virus
30641132	PVMoV - Positive control -	Pepper veinal mottle virus
30641133	PVMoV - Negative control -	Pepper veinal mottle virus
30641134	PVS - Reagent set -	Potato virus S
30641138	PVS - IDENTIKIT -	Potato virus S
30641139	PVS - TOTAL Kit -	Potato virus S
30641142	PVS - Positive control -	Potato virus S
30641143	PVS - Negative control -	Potato virus S
30641144	PVS - Reagent set -	Potato virus S
30641148	PVS - IDENTIKIT -	Potato virus S
30641149	PVS - TOTAL Kit -	Potato virus S
30641152	PVS - Positive control -	Potato virus S
30641153	PVS - Negative control -	Potato virus S
30641154	PVS - SPOTCHECK LF -	Potato virus S
30641156	PVT - Reagent set -	Potato virus T
30641160	PVT - IDENTIKIT -	Potato virus T
30641161	PVT - TOTAL Kit -	Potato virus T
30641164	PVT - Positive control -	Potato virus T
30641165	PVT - Negative control -	Potato virus T
30641166	PVV - Reagent set -	Potato virus V
30641170	PVV - IDENTIKIT -	Potato virus V
30641171	PVV - TOTAL Kit -	Potato virus V
30641174	PVV - Positive control -	Potato virus V
30641175	PVV - Negative control -	Potato virus V
30641176	PVV - SPOTCHECK LF -	Potato virus V

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641178	PVX - Reagent set -	Potato virus X - DAS
30641182	PVX - IDENTIKIT -	Potato virus X - DAS
30641183	PVX - TOTAL Kit -	Potato virus X - DAS
30641186	PVX - Positive control -	Potato virus X - DAS
30641187	PVX - Negative control -	Potato virus X - DAS
30641188	PVX - Reagent set -	Potato virus X - TAS
30641191	PVX - IDENTIKIT -	Potato virus X - TAS
30641192	PVX - TOTAL Kit -	Potato virus X - TAS
30641195	PVX - Positive control -	Potato virus X - TAS
30641196	PVX - Negative control -	Potato virus X - TAS
30641197	PVX - Reagent set -	Potato virus X
30641201	PVX - IDENTIKIT -	Potato virus X
30641202	PVX - TOTAL Kit -	Potato virus X
30641205	PVX - Positive control -	Potato virus X
30641206	PVX - Negative control -	Potato virus X
30641207	PVX - SPOTCHECK LF -	Potato virus X
30641209	PVYn - Reagent set -	Potato virus Yn - DAS
30641213	PVYn - IDENTIKIT -	Potato virus Yn - DAS
30641214	PVYn - TOTAL Kit -	Potato virus Yn - DAS
30641217	PVYn - Positive control -	Potato virus Yn - DAS
30641218	PVYn - Negative control	Potato virus Yn - DAS
30641219	PVYn - Reagent set -	Potato virus Yn - TAS
30641222	PVYn - IDENTIKIT -	Potato virus Yn - TAS
30641223	PVYn - TOTAL Kit -	Potato virus Yn - TAS
30641226	PVYn - Positive control -	Potato virus Yn - TAS
30641227	PVYn - Negative control -	Potato virus Yn - TAS
30641228	PVYo - Reagent set -	Potato virus Yo
30641231	PVYo - IDENTIKIT -	Potato virus Yo
30641232	PVYo - TOTAL Kit -	Potato virus Yo
30641235	PVYo - Positive control -	Potato virus Yo
30641236	PVYo - Negative control -	Potato virus Yo
30641237	PVYo/c - Reagent set -	Potato virus Yo/c
30641241	PVYo/c - IDENTIKIT -	Potato virus Yo/c
30641242	PVYo/c - TOTAL Kit -	Potato virus Yo/c
30641245	PVYo/c - Positive control -	Potato virus Yo/c
30641246	PVYo/c - Negative control -	Potato virus Yo/c

SKU	Product Description	Pathogen
30641247	PVYo/c/n - Reagent set -	Potato virus Yo/c/n - DAS
30641251	PVYo/c/n - IDENTIKIT -	Potato virus Yo/c/n - DAS
30641252	PVYo/c/n - TOTAL Kit -	Potato virus Yo/c/n - DAS
30641255	PVYo/c/n - Positive control -	Potato virus Yo/c/n - DAS
30641256	PVYo/c/n - Negative control -	Potato virus Yo/c/n - DAS
30641257	PVYo/c/n - Reagent set -	Potato virus Yo/c/n - TAS
30641260	PVYo/c/n - IDENTIKIT -	Potato virus Yo/c/n - TAS
30641261	PVYo/c/n - TOTAL Kit -	Potato virus Yo/c/n - TAS
30641264	PVYo/c/n - Positive control -	Potato virus Yo/c/n - TAS
30641265	PVYo/c/n - Negative control -	Potato virus Yo/c/n - TAS
30641266	PVY - Reagent set -	Potato virus Y
30641270	PVY - IDENTIKIT -	Potato virus Y
30641271	PVY - TOTAL Kit -	Potato virus Y
30641274	PVY - Positive control -	Potato virus Y
30641275	PVY - Negative control -	Potato virus Y
30641276	PVY - Reagent set -	Potato virus Y Analyser
30641279	PVY - IDENTIKIT -	Potato virus Y Analyser
30641280	PVY - TOTAL Kit -	Potato virus Y Analyser
30641283	PVY - Positive control -	Potato virus Y Analyser
30641284	PVY - Negative control -	Potato virus Y Analyser
30641285	PVY - SPOTCHECK LF -	Potato virus Y Analyser
30641287	PYDV - Reagent set -	Potato yellow dwarf virus
30641290	PYDV - IDENTIKIT -	Potato yellow dwarf virus
30641291	PYDV - TOTAL Kit -	Potato yellow dwarf virus
30641294	PYDV - Positive control -	Potato yellow dwarf virus
30641295	PYDV - Negative control -	Potato yellow dwarf virus
30641296	PZSV - Reagent set -	Pelargonium zonate spot virus
30641299	PZSV - IDENTIKIT -	Pelargonium zonate spot virus
30641300	PZSV - TOTAL Kit -	Pelargonium zonate spot virus
30641303	PZSV - Positive control -	Pelargonium zonate spot virus
30641304	PZSV - Negative control -	Pelargonium zonate spot virus
30641305	RaMV - Reagent set -	Radish mosaic virus
30641308	RaMV - IDENTIKIT -	Radish mosaic virus
30641309	RaMV - TOTAL Kit -	Radish mosaic virus
30641312	RaMV - Positive control -	Radish mosaic virus
30641313	RaMV - Negative control -	Radish mosaic virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641314	RBDV - Reagent set -	Raspberry bushy dwarf virus
30641317	RBDV - IDENTIKIT -	Raspberry bushy dwarf virus
30641318	RBDV - TOTAL Kit -	Raspberry bushy dwarf virus
30641321	RBDV - Positive control -	Raspberry bushy dwarf virus
30641322	RBDV - Negative control -	Raspberry bushy dwarf virus
30641323	RLSV - Reagent set -	Red lasoda virus (Potato latent virus)
30641326	RLSV - IDENTIKIT -	Red lasoda virus (Potato latent virus)
30641327	RLSV - TOTAL Kit -	Red lasoda virus (Potato latent virus)
30641330	RLSV - Positive control -	Red lasoda virus (Potato latent virus)
30641331	RLSV - Negative control -	Red lasoda virus (Potato latent virus)
30641332	RMV - Reagent set -	Ribgrass mosaic virus
30641335	RMV - IDENTIKIT -	Ribgrass mosaic virus
30641336	RMV - TOTAL Kit -	Ribgrass mosaic virus
30641339	RMV - Positive control -	Ribgrass mosaic virus
30641340	RMV - Negative control -	Ribgrass mosaic virus
30641341	RpRSV - Reagent set -	Raspberry ringspot virus
30641344	RpRSV - IDENTIKIT -	Raspberry ringspot virus
30641345	RpRSV - TOTAL Kit -	Raspberry ringspot virus
30641348	RpRSV - Positive control -	Raspberry ringspot virus
30641349	RpRSV - Negative control -	Raspberry ringspot virus
30641350	RSV - Reagent set -	Rice stripe virus
30641353	RSV - IDENTIKIT -	Rice stripe virus
30641354	RSV - TOTAL Kit -	Rice stripe virus
30641357	RSV - Positive control -	Rice stripe virus
30641358	RSV - Negative control -	Rice stripe virus
30641359	RYMV - Reagent set -	Rice yellow mottle virus
30641362	RYMV - IDENTIKIT -	Rice yellow mottle virus
30641363	RYMV - TOTAL Kit -	Rice yellow mottle virus
30641366	RYMV - Positive control -	Rice yellow mottle virus
30641367	RYMV - Negative control -	Rice yellow mottle virus
30641368	SBMV - Reagent set -	Southern bean mosaic virus
30641371	SBMV - IDENTIKIT -	Southern bean mosaic virus
30641372	SBMV - TOTAL Kit -	Southern bean mosaic virus
30641375	SBMV - Positive control -	Southern bean mosaic virus
30641376	SBMV - Negative control -	Southern bean mosaic virus
30641377	SBWMV - Reagent set -	Soil-borne wheat mosaic virus
30641380	SBWMV - IDENTIKIT -	Soil-borne wheat mosaic virus
30641381	SBWMV - TOTAL Kit -	Soil-borne wheat mosaic virus
30641384	SBWMV - Positive control -	Soil-borne wheat mosaic virus
30641385	SBWMV - Negative control -	Soil-borne wheat mosaic virus

Call 1.888.BIOPRO to order.

SKU	Product Description	Pathogen
30641386	SCMV - Reagent set -	Sugarcane mosaic virus
30641389	SCMV - IDENTIKIT -	Sugarcane mosaic virus
30641390	SCMV - TOTAL Kit -	Sugarcane mosaic virus
30641393	SCMV - Positive control -	Sugarcane mosaic virus
30641394	SCMV - Negative control -	Sugarcane mosaic virus
30641395	SqLCV - Reagent set -	Squash leaf curl virus
30641398	SqLCV - IDENTIKIT -	Squash leaf curl virus
30641399	SqLCV - TOTAL Kit	Squash leaf curl virus
30641402	SqLCV - Positive control -	Squash leaf curl virus
30641403	SqLCV - Negative control -	Squash leaf curl virus
30641404	SLRSV - Reagent set -	Strawberry latent ringspot virus
30641407	SLRSV - IDENTIKIT -	Strawberry latent ringspot virus
30641408	SLRSV - TOTAL Kit -	Strawberry latent ringspot virus
30641411	SLRSV - Positive control -	Strawberry latent ringspot virus
30641412	SLRSV - Negative control -	Strawberry latent ringspot virus
30641413	SLV - Reagent set -	Shallot latent virus
30641416	SLV - IDENTIKIT -	Shallot latent virus
30641417	SLV - TOTAL Kit -	Shallot latent virus
30641420	SLV - Positive control -	Shallot latent virus
30641421	SLV - Negative control -	Shallot latent virus
30641422	SMV - Reagent set -	Soybean mosaic virus
30641425	SMV - IDENTIKIT -	Soybean mosaic virus
30641426	SMV - TOTAL Kit -	Soybean mosaic virus
30641429	SMV - Positive control -	Soybean mosaic virus
30641430	SMV - Negative control -	Soybean mosaic virus
30641431	SMYEV - Reagent set -	Strawberry mild yellow edge virus
30641434	SMYEV - IDENTIKIT -	Strawberry mild yellow edge virus
30641435	SMYEV - TOTAL Kit -	Strawberry mild yellow edge virus
30641438	SMYEV - Positive control -	Strawberry mild yellow edge virus
30641439	SMYEV - Negative control -	Strawberry mild yellow edge virus
30641440	SoMV - Reagent set -	Sowbane mosaic virus
30641443	SoMV - IDENTIKIT -	Sowbane mosaic virus
30641444	SoMV - TOTAL Kit -	Sowbane mosaic virus
30641447	SoMV - Positive control -	Sowbane mosaic virus
30641448	SoMV - Negative control -	Sowbane mosaic virus
30641449	SPCSV - Reagent set -	Sweet potato chlorotic stunt virus
30641452	SPCSV - IDENTIKIT -	Sweet potato chlorotic stunt virus
30641453	SPCSV - TOTAL Kit -	Sweet potato chlorotic stunt virus
30641456	SPCSV - Positive control -	Sweet potato chlorotic stunt virus
30641457	SPCSV - Negative control -	Sweet potato chlorotic stunt virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641458	SPFMV - Reagent set -	Sweet potato feathery mottle virus
30641461	SPFMV - IDENTIKIT -	Sweet potato feathery mottle virus
30641462	SPFMV - TOTAL Kit -	Sweet potato feathery mottle virus
30641465	SPFMV - Positive control -	Sweet potato feathery mottle virus
30641466	SPFMV - Negative control -	Sweet potato feathery mottle virus
30641467	SqMV - Reagent set -	Squash mosaic virus
30641470	SqMV - IDENTIKIT -	Squash mosaic virus
30641471	SqMV - TOTAL Kit -	Squash mosaic virus
30641474	SqMV - Positive control -	Squash mosaic virus
30641475	SqMV - Negative control -	Squash mosaic virus
30641476	SrMV - Reagent set -	Sorghum mosaic virus
30641479	SrMV - IDENTIKIT -	Sorghum mosaic virus
30641480	SrMV - TOTAL Kit -	Sorghum mosaic virus
30641483	SrMV - Positive control -	Sorghum mosaic virus
30641484	SrMV - Negative control -	Sorghum mosaic virus
30641485	SYSV - Reagent set -	Shallot yellow stripe virus
30641488	SYSV - IDENTIKIT -	Shallot yellow stripe virus
30641489	SYSV - TOTAL Kit -	Shallot yellow stripe virus
30641492	SYSV - Positive control -	Shallot yellow stripe virus
30641493	SYSV - Negative control -	Shallot yellow stripe virus
30641494	TAV - Reagent set -	Tomato aspermy virus
30641497	TAV - IDENTIKIT -	Tomato aspermy virus
30641498	TAV - TOTAL Kit -	Tomato aspermy virus
30641501	TAV - Positive control -	Tomato aspermy virus
30641502	TAV - Negative control -	Tomato aspermy virus
30641503	TBRV - Reagent set -	Tomato blackring virus
30641506	TBRV - IDENTIKIT -	Tomato blackring virus
30641507	TBRV - TOTAL Kit -	Tomato blackring virus
30641510	TBRV - Positive control -	Tomato blackring virus
30641511	TBRV - Negative control -	Tomato blackring virus
30641512	TBSV - Reagent set -	Tomato bushy stunt virus
30641515	TBSV - IDENTIKIT -	Tomato bushy stunt virus
30641516	TBSV - TOTAL Kit -	Tomato bushy stunt virus
30641519	TBSV - Positive control -	Tomato bushy stunt virus
30641520	TBSV - Negative control -	Tomato bushy stunt virus
30641521	ToCSV - Reagent set -	Tomato chlorotic spot virus
30641524	ToCSV - IDENTIKIT -	Tomato chlorotic spot virus
30641525	ToCSV - TOTAL Kit -	Tomato chlorotic spot virus
30641528	ToCSV - Positive control -	Tomato chlorotic spot virus
30641529	ToCSV - Negative control -	Tomato chlorotic spot virus

Call 1.888.BIOPRO to order.

SKU	Product Description	Pathogen
30641530	TEV - Reagent set -	Tobacco etch virus
30641533	TEV - IDENTIKIT -	Tobacco etch virus
30641534	TEV - TOTAL Kit -	Tobacco etch virus
30641537	TEV - Positive control -	Tobacco etch virus
30641538	TEV - Negative control -	Tobacco etch virus
30641539	TMV - Reagent set -	Tobacco mosaic virus
30641542	TMV - IDENTIKIT -	Tobacco mosaic virus
30641543	TMV - TOTAL Kit -	Tobacco mosaic virus
30641546	TMV - Positive control -	Tobacco mosaic virus
30641547	TMV - Negative control -	Tobacco mosaic virus
30641548	TNV - Reagent set -	Tobacco necrosis virus
30641551	TNV - IDENTIKIT -	Tobacco necrosis virus
30641552	TNV - TOTAL Kit -	Tobacco necrosis virus
30641555	TNV - Positive control -	Tobacco necrosis virus
30641556	TNV - Negative control -	Tobacco necrosis virus
30641557	ToMV - Reagent set -	Tomato mosaic virus
30641560	ToMV - IDENTIKIT -	Tomato mosaic virus
30641561	ToMV - TOTAL Kit -	Tomato mosaic virus
30641564	ToMV - Positive control -	Tomato mosaic virus
30641565	ToMV - Negative control -	Tomato mosaic virus
30641566	ToMV - SPOTCHECK LF -	Tomato mosaic virus
30641568	ToRSV - Reagent set -	Tomato ringspot virus
30641571	ToRSV - IDENTIKIT -	Tomato ringspot virus
30641572	ToRSV - TOTAL Kit -	Tomato ringspot virus
30641575	ToRSV - Positive control -	Tomato ringspot virus
30641576	ToRSV - Negative control -	Tomato ringspot virus
30641577	TRSV - Reagent set -	Tobacco ringspot virus - Calico strain
30641580	TRSV - IDENTIKIT -	Tobacco ringspot virus - Calico strain
30641581	TRSV - TOTAL Kit -	Tobacco ringspot virus - Calico strain
30641584	TRSV - Positive control -	Tobacco ringspot virus - Calico strain
30641585	TRSV - Negative control -	Tobacco ringspot virus - Calico strain
30641586	TRV - Reagent set -	Tobacco rattle virus
30641589	TRV - IDENTIKIT -	Tobacco rattle virus
30641590	TRV - TOTAL Kit -	Tobacco rattle virus
30641593	TRV - Positive control -	Tobacco rattle virus
30641594	TRV - Negative control -	Tobacco rattle virus

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641595	TSV - Reagent set -	Tobacco streak virus
30641598	TSV - IDENTIKIT -	Tobacco streak virus
30641599	TSV - TOTAL Kit -	Tobacco streak virus
30641602	TSV - Positive control -	Tobacco streak virus
30641603	TSV - Negative control -	Tobacco streak virus
30641604	TSWV - Reagent set -	Tomato spotted wilt virus
30641607	TSWV - IDENTIKIT -	Tomato spotted wilt virus
30641608	TSWV - TOTAL Kit -	Tomato spotted wilt virus
30641611	TSWV - Positive control -	Tomato spotted wilt virus
30641612	TSWV - Negative control -	Tomato spotted wilt virus
30641613	TSWV - SPOTCHECK LF -	Tomato spotted wilt virus
30641615	TuMV - Reagent set -	Turnip mosaic virus
30641618	TuMV - IDENTIKIT -	Turnip mosaic virus
30641619	TuMV - TOTAL Kit -	Turnip mosaic virus
30641622	TuMV - Positive control -	Turnip mosaic virus
30641623	TuMV - Negative control -	Turnip mosaic virus
30641624	TuYMV - Reagent set -	Turnip yellow mosaic virus
30641627	TuYMV - IDENTIKIT -	Turnip yellow mosaic virus
30641628	TuYMV - TOTAL Kit -	Turnip yellow mosaic virus
30641631	TuYMV - Positive control -	Turnip yellow mosaic virus
30641632	TuYMV - Negative control -	Turnip yellow mosaic virus
30641633	TVBMV - Reagent set -	Tobacco vein banding mosaic virus
30641636	TVBMV - IDENTIKIT -	Tobacco vein banding mosaic virus
30641637	TVBMV - TOTAL Kit -	Tobacco vein banding mosaic virus
30641640	TVBMV - Positive control -	Tobacco vein banding mosaic virus
30641641	TVBMV - Negative control -	Tobacco vein banding mosaic virus
30641642	TYLCV - Reagent set -	Tomato yellow leaf curl virus
30641645	TYLCV - IDENTIKIT -	Tomato yellow leaf curl virus
30641646	TYLCV - TOTAL Kit -	Tomato yellow leaf curl virus
30641649	TYLCV - Positive control -	Tomato yellow leaf curl virus
30641650	TYLCV - Negative control -	Tomato yellow leaf curl virus
30641651	WDV - Reagent set -	Wheat dwarf virus
30641654	WDV - IDENTIKIT -	Wheat dwarf virus
30641655	WDV - TOTAL Kit -	Wheat dwarf virus
30641658	WDV - Positive control -	Wheat dwarf virus
30641659	WDV - Negative control -	Wheat dwarf virus

SKU	Product Description	Pathogen
30641660	WMV-2 - Reagent set -	Watermelon mosaic virus - 2
30641663	WMV-2 - IDENTIKIT -	Watermelon mosaic virus - 2
30641664	WMV-2 - TOTAL Kit -	Watermelon mosaic virus - 2
30641667	WMV-2 - Positive control -	Watermelon mosaic virus - 2
30641668	WMV-2 - Negative control -	Watermelon mosaic virus - 2
30641669	WSMoV - Reagent set -	Watermelon silver mottle virus
30641672	WSMoV - IDENTIKIT -	Watermelon silver mottle virus
30641673	WSMoV - TOTAL Kit -	Watermelon silver mottle virus
30641676	WSMoV - Positive control -	Watermelon silver mottle virus
30641677	WSMoV - Negative control -	Watermelon silver mottle virus
30641678	WSMV - Reagent set -	Wheat streak mosaic virus
30641681	WSMV - IDENTIKIT -	Wheat streak mosaic virus
30641682	WSMV - TOTAL Kit -	Wheat streak mosaic virus
30641685	WSMV - Positive control -	Wheat streak mosaic virus
30641686	WSMV - Negative control -	Wheat streak mosaic virus
30641687	WSSMV - Reagent set -	Wheat spindle streak mosaic virus
30641690	WSSMV - IDENTIKIT -	Wheat spindle streak mosaic virus
30641691	WSSMV - TOTAL Kit -	Wheat spindle streak mosaic virus
30641694	WSSMV - Positive control -	Wheat spindle streak mosaic virus
30641695	WSSMV - Negative control -	Wheat spindle streak mosaic virus
30641696	ZGMMV - Reagent set -	Zucchini green mottle mosaic virus
30641699	ZGMMV - IDENTIKIT -	Zucchini green mottle mosaic virus
30641700	ZGMMV - TOTAL Kit -	Zucchini green mottle mosaic virus
30641703	ZGMMV - Positive control -	Zucchini green mottle mosaic virus
30641704	ZGMMV - Negative control -	Zucchini green mottle mosaic virus
30641705	ZYMV - Reagent set -	Zucchini yellow mosaic virus
30641708	ZYMV - IDENTIKIT -	Zucchini yellow mosaic virus
30641709	ZYMV - TOTAL Kit -	Zucchini yellow mosaic virus
30641712	ZYMV - Positive control -	Zucchini yellow mosaic virus
30641713	ZYMV - Negative control -	Zucchini yellow mosaic virus
30641714	B. andropogonis - Reagent set -	Burkholderia andropogonis
30641717	B. andropogonis - IDENTIKIT -	Burkholderia andropogonis
30641718	B. andropogonis - TOTAL Kit -	Burkholderia andropogonis
30641721	B. andropogonis - Positive control -	Burkholderia andropogonis
30641722	B. andropogonis - Negative control -	Burkholderia andropogonis
30641723	B. andropogonis - EXPRESS Kit -	Burkholderia andropogonis
30641725	B. andropogonis - Fluoriscan-IF Kit -	Burkholderia andropogonis

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641728	B. g. pv alliicola - Reagent set -	Burkholderia gladioli pv alliicola
30641731	B. g. pv alliicola - IDENTIKIT -	Burkholderia gladioli pv alliicola
30641732	B. g. pv alliicola - TOTAL Kit -	Burkholderia gladioli pv alliicola
30641735	B. g. pv alliicola - Positive control -	Burkholderia gladioli pv alliicola
30641736	B. g. pv alliicola - Negative control -	Burkholderia gladioli pv alliicola
30641737	B. g. pv alliicola - EXPRESS Kit -	Burkholderia gladioli pv alliicola
30641739	B. g. pv alliicola - Fluoriscan-IF Kit -	Burkholderia gladioli pv alliicola
30641742	C. m. subsp. insidiosum - Reagent set -	Clavibacter michiganensis subsp. insidiosum
30641745	C. m. subsp. insidiosum - IDENTIKIT -	Clavibacter michiganensis subsp. insidiosum
30641746	C. m. subsp. insidiosum - TOTAL Kit -	Clavibacter michiganensis subsp. insidiosum
30641749	C. m. subsp. insidiosum - Positive control -	Clavibacter michiganensis subsp. insidiosum
30641750	C. m. subsp. insidiosum - Negative control -	Clavibacter michiganensis subsp. insidiosum
30641751	C. m. subsp. insidiosum - EXPRESS Kit -	Clavibacter michiganensis subsp. insidiosum
30641753	C. m. subsp. insidiosum - Fluoriscan-IF Kit -	Clavibacter michiganensis subsp. insidiosum
30641756	C. m. subsp. michiganensis - Reagent set -	Clavibacter michiganensis subsp. michiganensis
30641759	C. m. subsp. michiganensis - IDENTIKIT -	Clavibacter michiganensis subsp. michiganensis
30641760	C. m. subsp. michiganensis - TOTAL Kit -	Clavibacter michiganensis subsp. michiganensis
30641763	C. m. subsp. michiganensis - Positive control -	Clavibacter michiganensis subsp. michiganensis
30641764	C. m. subsp. michiganensis - Negative control -	Clavibacter michiganensis subsp. michiganensis
30641765	C. m. subsp. michiganensis - EXPRESS Kit -	Clavibacter michiganensis subsp. michiganensis
30641767	C. m. subsp. michiganensis - Fluoriscan-IF Kit -	Clavibacter michiganensis subsp. michiganensis
30641770	C. m. subsp. nebraskensis - Reagent set -	Clavibacter michiganensis subsp. nebraskensis
30641773	C. m. subsp. nebraskensis - IDENTIKIT -	Clavibacter michiganensis subsp. nebraskensis
30641774	C. m. subsp. nebraskensis - TOTAL Kit -	Clavibacter michiganensis subsp. nebraskensis
30641777	C. m. subsp. nebraskensis - Positive control -	Clavibacter michiganensis subsp. nebraskensis
30641778	C. m. subsp. nebraskensis - Negative control -	Clavibacter michiganensis subsp. nebraskensis
30641779	C. m. subsp. sepidonicus - Reagent set -	Clavibacter michiganensis subsp. sepidonicus
30641782	C. m. subsp. sepidonicus - IDENTIKIT -	Clavibacter michiganensis subsp. sepidonicus
30641783	C. m. subsp. sepidonicus - TOTAL Kit -	Clavibacter michiganensis subsp. sepidonicus
30641786	C. m. subsp. sepidonicus - Positive control -	Clavibacter michiganensis subsp. sepidonicus
30641787	C. m. subsp. sepidonicus - Negative control -	Clavibacter michiganensis subsp. sepidonicus
30641788	C. m. subsp. sepidonicus - EXPRESS Kit -	Clavibacter michiganensis subsp. sepidonicus
30641790	C. m. subsp. sepidonicus - Fluoriscan-IF Kit -	Clavibacter michiganensis subsp. sepidonicus
30641793	C. f. subsp. oortii - Reagent set -	Curtobacterium flaccumfaciens subsp. oortii
30641796	C. f. subsp. oortii - IDENTIKIT -	Curtobacterium flaccumfaciens subsp. oortii
30641797	C. f. subsp. oortii - TOTAL Kit -	Curtobacterium flaccumfaciens subsp. oortii
30641800	C. f. subsp. oortii - Positive control -	Curtobacterium flaccumfaciens subsp. oortii
30641801	C. f. subsp. oortii - Negative control -	Curtobacterium flaccumfaciens subsp. oortii
30641802	C. f. subsp. oortii - EXPRESS Kit -	Curtobacterium flaccumfaciens subsp. oortii
30641804	C. f. subsp. oortii - Fluoriscan-IF Kit -	Curtobacterium flaccumfaciens subsp. oortii

Call 1.888.BIOPPLUS to order.

SKU	Product Description	Pathogen
30641807	C. f. subsp. flaccumfaciens - Reagent set -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641810	C. f. subsp. flaccumfaciens - IDENTIKIT -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641811	C. f. subsp. flaccumfaciens - TOTAL Kit -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641814	C. f. subsp. flaccumfaciens - Positive control -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641815	C. f. subsp. flaccumfaciens - Negative control -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641816	C. f. subsp. flaccumfaciens - EXPRESS Kit -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641818	C. f. subsp. flaccumfaciens - Fluoriscan-IF Kit -	<i>Curtobacterium flaccumfaciens</i> subsp. <i>flaccumfaciens</i>
30641821	E. amylovora - Reagent set -	<i>Erwinia amylovora</i>
30641824	E. amylovora - IDENTIKIT -	<i>Erwinia amylovora</i>
30641825	E. amylovora - TOTAL Kit -	<i>Erwinia amylovora</i>
30641828	E. amylovora - Positive control -	<i>Erwinia amylovora</i>
30641829	E. amylovora - Negative control -	<i>Erwinia amylovora</i>
30641830	E. amylovora - EXPRESS Kit -	<i>Erwinia amylovora</i>
30641832	E. amylovora - Fluoriscan-IF Kit -	<i>Erwinia amylovora</i>
30641835	E. amylovora - Reagent set -	<i>Erwinia amylovora</i>
30641838	E. amylovora - IDENTIKIT -	<i>Erwinia amylovora</i>
30641841	E. amylovora - TOTAL Kit -	<i>Erwinia amylovora</i>
30641842	E. amylovora - Positive control -	<i>Erwinia amylovora</i>
30641843	E. amylovora - Negative control -	<i>Erwinia amylovora</i>
30641844	E. amylovora - SPOTCHECK LF -	<i>Erwinia amylovora</i>
30641846	E. c. subsp. atroseptica - Reagent set -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641849	E. c. subsp. atroseptica - IDENTIKIT -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641850	E. c. subsp. atroseptica - TOTAL Kit -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641853	E. c. subsp. atroseptica - Positive control -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641854	E. c. subsp. atroseptica - Negative control -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641855	E. c. subsp. atroseptica - EXPRESS Kit -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641857	E. c. subsp. atroseptica - Fluoriscan-IF Kit -	<i>Erwinia carotovora</i> subsp. <i>atroseptica</i>
30641860	E. chrysanthemi - Reagent set -	<i>Erwinia chrysanthemi</i>
30641863	E. chrysanthemi - IDENTIKIT -	<i>Erwinia chrysanthemi</i>
30641864	E. chrysanthemi - TOTAL Kit -	<i>Erwinia chrysanthemi</i>
30641867	E. chrysanthemi - Positive control -	<i>Erwinia chrysanthemi</i>
30641868	E. chrysanthemi - Negative control -	<i>Erwinia chrysanthemi</i>
30641869	E. chrysanthemi - EXPRESS Kit -	<i>Erwinia chrysanthemi</i>
30641871	E. chrysanthemi - Fluoriscan-IF Kit -	<i>Erwinia chrysanthemi</i>
30641874	E. stewartii - Reagent set -	<i>Erwinia stewartii</i>
30641877	E. stewartii - IDENTIKIT -	<i>Erwinia stewartii</i>
30641878	E. stewartii - TOTAL Kit -	<i>Erwinia stewartii</i>
30641881	E. stewartii - Positive control -	<i>Erwinia stewartii</i>
30641882	E. stewartii - Negative control -	<i>Erwinia stewartii</i>
30641883	E. stewartii - EXPRESS Kit -	<i>Erwinia stewartii</i>
30641885	E. stewartii - Fluoriscan-IF Kit -	<i>Erwinia stewartii</i>

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30641888	F. doree - Reagent set -	Flavescence doree
30641891	F. doree - IDENTIKIT -	Flavescence doree
30641892	F. doree - TOTAL Kit -	Flavescence doree
30641895	F. doree - Positive control -	Flavescence doree
30641896	F. doree - Negative control -	Flavescence doree
30641897	P. s. pv glycinea - Reagent set -	Pseudomonas syringae pv glycinea
30641900	P. s. pv glycinea - IDENTIKIT -	Pseudomonas syringae pv glycinea
30641901	P. s. pv glycinea - TOTAL Kit -	Pseudomonas syringae pv glycinea
30641904	P. s. pv glycinea - Positive control -	Pseudomonas syringae pv glycinea
30641905	P. s. pv glycinea - Negative control -	Pseudomonas syringae pv glycinea
30641906	P. s. pv glycinea - EXPRESS Kit -	Pseudomonas syringae pv glycinea
30641908	P. s. pv glycinea - Fluoriscan-IF Kit -	Pseudomonas syringae pv glycinea
30641925	P. s. pv phaseolicola - Reagent set -	Pseudomonas syringae pv phaseolicola
30641928	P. s. pv phaseolicola - IDENTIKIT -	Pseudomonas syringae pv phaseolicola
30641929	P. s. pv phaseolicola - TOTAL Kit -	Pseudomonas syringae pv phaseolicola
30641932	P. s. pv phaseolicola - Positive control -	Pseudomonas syringae pv phaseolicola
30641933	P. s. pv phaseolicola - Negative control -	Pseudomonas syringae pv phaseolicola
30641934	P. s. pv phaseolicola - EXPRESS Kit -	Pseudomonas syringae pv phaseolicola
30641936	P. s. pv phaseolicola - Fluoriscan-IF Kit -	Pseudomonas syringae pv phaseolicola
30641925	P. s. pv phaseolicola - Reagent set -	Pseudomonas syringae pv phaseolicola
30641928	P. s. pv phaseolicola - IDENTIKIT -	Pseudomonas syringae pv phaseolicola
30641929	P. s. pv phaseolicola - TOTAL Kit -	Pseudomonas syringae pv phaseolicola
30641932	P. s. pv phaseolicola - Positive control -	Pseudomonas syringae pv phaseolicola
30641933	P. s. pv phaseolicola - Negative control -	Pseudomonas syringae pv phaseolicola
30641934	P. s. pv phaseolicola - EXPRESS Kit -	Pseudomonas syringae pv phaseolicola
30641936	P. s. pv phaseolicola - Fluoriscan-IF Kit -	Pseudomonas syringae pv phaseolicola
30641939	P. s. pv tomato - Reagent set -	Pseudomonas syringae pv tomato
30641942	P. s. pv tomato - IDENTIKIT -	Pseudomonas syringae pv tomato
30641943	P. s. pv tomato - TOTAL Kit -	Pseudomonas syringae pv tomato
30641946	P. s. pv tomato - Positive control -	Pseudomonas syringae pv tomato
30641947	P. s. pv tomato - Negative control -	Pseudomonas syringae pv tomato
30641948	P. s. pv tomato - EXPRESS Kit -	Pseudomonas syringae pv tomato
30641950	P. s. pv tomato - Fluoriscan-IF Kit -	Pseudomonas syringae pv tomato
30641953	P. s. pv pisi - Reagent set -	Pseudomonas syringae pv pisi
30641956	P. s. pv pisi - IDENTIKIT -	Pseudomonas syringae pv pisi
30641957	P. s. pv pisi - TOTAL Kit -	Pseudomonas syringae pv pisi
30641960	P. s. pv pisi - Positive control -	Pseudomonas syringae pv pisi
30641961	P. s. pv pisi - Negative control -	Pseudomonas syringae pv pisi
30641962	P. s. pv pisi - EXPRESS Kit -	Pseudomonas syringae pv pisi
30641964	P. s. pv pisi - Fluoriscan-IF Kit -	Pseudomonas syringae pv pisi

SKU	Product Description	Pathogen
30641967	P. s. pv syringae - Reagent set -	Pseudomonas syringae pv syringae
30641970	P. s. pv syringae - IDENTIKIT -	Pseudomonas syringae pv syringae
30641971	P. s. pv syringae - TOTAL Kit -	Pseudomonas syringae pv syringae
30641974	P. s. pv syringae - Positive control -	Pseudomonas syringae pv syringae
30641975	P. s. pv syringae - Negative control - 2ml	Pseudomonas syringae pv syringae
30641976	P. s. pv syringae - EXPRESS Kit - 50 tests	Pseudomonas syringae pv syringae
30641978	P. s. pv syringae - Fluoriscan-IF Kit - 250 tests	Pseudomonas syringae pv syringae
30641981	R. solonacearum - Reagent set - 250 units	Ralstonia solonacearum
30641984	R. solonacearum - IDENTIKIT -	Ralstonia solonacearum
30641985	R. solonacearum - TOTAL Kit -	Ralstonia solonacearum
30641988	R. solonacearum - Positive control -	Ralstonia solonacearum
30641989	R. solonacearum - Negative control -	Ralstonia solonacearum
30641990	R. solonacearum - EXPRESS Kit -	Ralstonia solonacearum
30641992	R. solonacearum - Fluoriscan-IF Kit -	Ralstonia solonacearum
30641997	X. c. pv begoniae - Reagent set -	Xanthomonas campestris pv begoniae
30642000	X. c. pv begoniae - IDENTIKIT -	Xanthomonas campestris pv begoniae
30642001	X. c. pv begoniae - TOTAL Kit -	Xanthomonas campestris pv begoniae
30642004	X. c. pv begoniae - Positive control -	Xanthomonas campestris pv begoniae
30642005	X. c. pv begoniae - Negative control -	Xanthomonas campestris pv begoniae
30642006	X. c. pv carotae - Reagent set -	Xanthomonas campestris pv carotae
30642009	X. c. pv carotae - IDENTIKIT -	Xanthomonas campestris pv carotae
30642010	X. c. pv carotae - TOTAL Kit -	Xanthomonas campestris pv carotae
30642013	X. c. pv carotae - Positive control -	Xanthomonas campestris pv carotae
30642014	X. c. pv carotae - Negative control -	Xanthomonas campestris pv carotae
30642015	X. c. pv carotae - EXPRESS Kit -	Xanthomonas campestris pv carotae
30642017	X. c. pv carotae - Fluoriscan-IF Kit -	Xanthomonas campestris pv carotae
30642020	X. c. pv campestris - Reagent set -	Xanthomonas campestris pv campestris
30642023	X. c. pv campestris - IDENTIKIT -	Xanthomonas campestris pv campestris
30642024	X. c. pv campestris - TOTAL Kit -	Xanthomonas campestris pv campestris
30642027	X. c. pv campestris - Positive control -	Xanthomonas campestris pv campestris
30642028	X. c. pv campestris - Negative control -	Xanthomonas campestris pv campestris
30642029	X. c. pv campestris - EXPRESS Kit -	Xanthomonas campestris pv campestris
30642031	X. c. pv campestris - Fluoriscan-IF Kit -	Xanthomonas campestris pv campestris
30642034	X. c. pv dieffenbachiae - Reagent set -	Xanthomonas campestris pv dieffenbachiae
30642037	X. c. pv dieffenbachiae - IDENTIKIT -	Xanthomonas campestris pv dieffenbachiae
30642038	X. c. pv dieffenbachiae - TOTAL Kit - 250 units	Xanthomonas campestris pv dieffenbachiae
30642041	X. c. pv dieffenbachiae - Positive control - 2ml	Xanthomonas campestris pv dieffenbachiae
30642042	X. c. pv dieffenbachiae - Negative control - 2ml	Xanthomonas campestris pv dieffenbachiae
30642043	X. c. pv dieffenbachiae - EXPRESS Kit - 50 tests	Xanthomonas campestris pv dieffenbachiae
30642045	X. c. pv dieffenbachiae - Fluoriscan-IF Kit - 250 tests	Xanthomonas campestris pv dieffenbachiae

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30642048	X. c. pv oryzicola - Reagent set -	Xanthomonas campestris pv oryzicola
30642051	X. c. pv oryzicola - IDENTIKIT -	Xanthomonas campestris pv oryzicola
30642052	X. c. pv oryzicola - TOTAL Kit -	Xanthomonas campestris pv oryzicola
30642055	X. c. pv oryzicola - Positive control -	Xanthomonas campestris pv oryzicola
30642056	X. c. pv oryzicola - Negative control -	Xanthomonas campestris pv oryzicola
30642057	X. c. pv oryzicola - EXPRESS Kit -	Xanthomonas campestris pv oryzicola
30642059	X. c. pv oryzicola - Fluoriscan-IF Kit -	Xanthomonas campestris pv oryzicola
30642062	X. c. pv pelargonii - Reagent set -	Xanthomonas campestris pv pelargonii
30642065	X. c. pv pelargonii - IDENTIKIT -	Xanthomonas campestris pv pelargonii
30642066	X. c. pv pelargonii - TOTAL Kit -	Xanthomonas campestris pv pelargonii
30642069	X. c. pv pelargonii - Positive control -	Xanthomonas campestris pv pelargonii
30642070	X. c. pv pelargonii - Negative control -	Xanthomonas campestris pv pelargonii
30642071	X. c. pv pelargonii - EXPRESS Kit -	Xanthomonas campestris pv pelargonii
30642073	X. c. pv pelargonii - Fluoriscan-IF Kit -	Xanthomonas campestris pv pelargonii
30642078	X. c. pv phaseoli - Reagent set -	Xanthomonas campestris pv phaseoli
30642081	X. c. pv phaseoli - IDENTIKIT -	Xanthomonas campestris pv phaseoli
30642082	X. c. pv phaseoli - TOTAL Kit -	Xanthomonas campestris pv phaseoli
30642085	X. c. pv phaseoli - Positive control -	Xanthomonas campestris pv phaseoli
30642086	X. c. pv phaseoli - Negative control -	Xanthomonas campestris pv phaseoli
30642087	X. c. pv phaseoli - EXPRESS Kit -	Xanthomonas campestris pv phaseoli
30642089	X. c. pv phaseoli - Fluoriscan-IF Kit -	Xanthomonas campestris pv phaseoli
30642092	X. c. pv pruni - Reagent set -	Xanthomonas campestris pv pruni
30642095	X. c. pv pruni - IDENTIKIT -	Xanthomonas campestris pv pruni
30642096	X. c. pv pruni - TOTAL Kit -	Xanthomonas campestris pv pruni
30642099	X. c. pv pruni - Positive control -	Xanthomonas campestris pv pruni
30642100	X. c. pv pruni - Negative control -	Xanthomonas campestris pv pruni
30642101	X. c. pv pruni - EXPRESS Kit -	Xanthomonas campestris pv pruni
30642103	X. c. pv pruni - Fluoriscan-IF Kit -	Xanthomonas campestris pv pruni
30642106	X. c. pv vesicatoria - Reagent set -	Xanthomonas campestris pv vesicatoria
30642109	X. c. pv vesicatoria - IDENTIKIT -	Xanthomonas campestris pv vesicatoria
30642110	X. c. pv vesicatoria - TOTAL Kit -	Xanthomonas campestris pv vesicatoria
30642113	X. c. pv vesicatoria - Positive control -	Xanthomonas campestris pv vesicatoria
30642114	X. c. pv vesicatoria - Negative control -	Xanthomonas campestris pv vesicatoria
30642115	X. c. pv vesicatoria - EXPRESS Kit -	Xanthomonas campestris pv vesicatoria
30642117	X. c. pv vesicatoria - Fluoriscan-IF Kit -	Xanthomonas campestris pv vesicatoria

SKU	Product Description	Pathogen
30642120	X. fragariae - Reagent set -	Xanthomonas fragariae
30642123	X. fragariae - IDENTIKIT -	Xanthomonas fragariae
30642124	X. fragariae - TOTAL Kit -	Xanthomonas fragariae
30642127	X. fragariae - Positive control -	Xanthomonas fragariae
30642128	X. fragariae - Negative control -	Xanthomonas fragariae
30642129	X. fragariae - EXPRESS Kit -	Xanthomonas fragariae
30642131	X. fragariae - Fluoriscan-IF Kit -	Xanthomonas fragariae
30642134	B. allii - Reagent set -	Botrytis allii
30642137	B. allii - IDENTIKIT -	Botrytis allii
30642138	B. allii - TOTAL Kit -	Botrytis allii
30642141	B. allii - Positive control -	Botrytis allii
30642142	B. allii - Negative control -	Botrytis allii
30642143	B. cinerea - Reagent set -	Botrytis cinerea
30642146	B. cinerea - IDENTIKIT-Q -	Botrytis cinerea
30642147	B. cinerea - TOTAL Kit -	Botrytis cinerea
30642150	B. cinerea - Positive control -	Botrytis cinerea
30642151	B. cinerea - Negative control -	Botrytis cinerea
30642152	Botrytis spp. - ALERT LF -	Botrytis spp.
30642154	C. acutatum - Reagent set -	Colletotrichum acutatum
30642157	C. acutatum - IDENTIKIT -	Colletotrichum acutatum
30642158	C. acutatum - TOTAL Kit -	Colletotrichum acutatum
30642161	C. acutatum - Positive control -	Colletotrichum acutatum
30642162	C. acutatum - Negative control -	Colletotrichum acutatum
30642163	Colletotrichum spp. - Reagent set -	Colletotrichum spp.
30642166	Colletotrichum spp. - IDENTIKIT -	Colletotrichum spp.
30642167	Colletotrichum spp. - TOTAL Kit -	Colletotrichum spp.
30642170	Colletotrichum spp. - Positive control -	Colletotrichum spp.
30642171	Colletotrichum spp. - Negative control -	Colletotrichum spp.
30642172	Phytophthora spp. - Reagent set -	Phytophthora spp.
30642175	Phytophthora spp. - AGRISCREEN -	Phytophthora spp.
30642176	Phytophthora spp. - TOTAL Kit -	Phytophthora spp.
30642179	Phytophthora spp. - Positive control -	Phytophthora spp.
30642180	Phytophthora spp. - Negative control -	Phytophthora spp.
30642181	Phytophthora spp. - ALERT LF -	Phytophthora spp.

For available sizes & current prices, visit: www.plantmedia.com

SKU	Product Description	Pathogen
30642183	Pythium spp. - Reagent set -	Pythium spp.
30642186	Pythium spp. - AGRISCREEN -	Pythium spp.
30642187	Pythium spp. - TOTAL Kit -	Pythium spp.
30642190	Pythium spp. - Positive control -	Pythium spp.
30642191	Pythium spp. - Negative control -	Pythium spp.
30642192	Pythium spp. - ALERT LF -	Pythium spp.
30642194	Rhizoctonia spp. - Reagent set -	Rhizoctonia spp.
30642197	Rhizoctonia spp. - AGRISCREEN - 96 wells	Rhizoctonia spp.
30642198	Rhizoctonia spp. - TOTAL Kit -	Rhizoctonia spp.
30642201	Rhizoctonia spp. - Positive control -	Rhizoctonia spp.
30642202	Rhizoctonia spp. - Negative control -	Rhizoctonia spp.
30642203	Rhizoctonia spp. - ALERT LF -	Rhizoctonia spp.
30642205	S. tritici - Reagent set -	Septoria tritici
30642208	S. tritici - IDENTIKIT-Q -	Septoria tritici
30642209	S. tritici - TOTAL Kit -	Septoria tritici
30642212	S. tritici - Positive control -	Septoria tritici
30642213	S. tritici - Negative control -	Septoria tritici
30642214	S. nodorum / S. tritici - Reagent set -	Septoria nodorum / Septoria tritici
30642217	S. nodorum / S. tritici - IDENTIKIT-Q -	Septoria nodorum / Septoria tritici
30642218	S. nodorum / S. tritici - TOTAL Kit -	Septoria nodorum / Septoria tritici
30642221	S. nodorum / S. tritici - Positive control -	Septoria nodorum / Septoria tritici
30642222	S. nodorum / S. tritici - Negative control -	Septoria nodorum / Septoria tritici
30642223	Tapesia spp. - Reagent set -	Tapesia spp.
30642226	Tapesia spp. - IDENTIKIT-Q -	Tapesia spp.
30642227	Tapesia spp. - TOTAL Kit -	Tapesia spp.
30642230	Tapesia spp. - Positive control -	Tapesia spp.
30642231	Tapesia spp. - Negative control -	Tapesia spp.



Sunrise Science products now available at bio-world.com

Yeast Dropout Media


SKU #	Product Name	Product Group
30626100	CSM Powder, 10 grams	0-amino acid dropouts
30626101	CSM-His-Leu-Trp-Ura Powder, 10 grams	4-amino acid dropouts
30626102	CSM-Ura Powder, 10 grams	1-amino acid dropouts
30626103	CSM-Leu Powder, 10 grams	1-amino acid dropouts
30626104	CSM-His Powder, 10 grams	1-amino acid dropouts
30626105	CSM-Trp Powder, 10 grams	1-amino acid dropouts
30626106	CSM-His-Ura Powder, 10 grams	2-amino acid dropouts
30626107	CSM-Trp-Ura Powder, 10 grams	2-amino acid dropouts
30626108	CSM-His-Leu Powder, 10 grams	2-amino acid dropouts
30626109	CSM-Leu-Trp Powder, 10 grams	2-amino acid dropouts
30626110	CSM-His-Trp Powder, 10 grams	2-amino acid dropouts
30626111	CSM-His-Leu-Trp Powder, 10 grams	3-amino acid dropouts
30626112	CSM-His-Leu-Ura Powder, 10 grams	3-amino acid dropouts
30626113	CSM-His-Trp-Ura Powder, 10 grams	3-amino acid dropouts
30626114	CSM-Leu-Trp-Ura Powder, 10 grams	3-amino acid dropouts
30626115	CSM-Lys-Ura Powder, 10 grams	2-amino acid dropouts
30626116	CSM-Arg-Ura Powder, 10 grams	2-amino acid dropouts
30626117	CSM-Ade-His-Leu Powder, 10 grams	3-amino acid dropouts
30626118	CSM-Ade-Leu-Trp Powder, 10 grams	3-amino acid dropouts
30626119	CSM-Ade-His-Leu-Trp Powder, 10 grams	4-amino acid dropouts
30626120	CSM-Ade-Trp Powder, 10 grams	2-amino acid dropouts
30626121	CSM-Ade Powder, 10 grams	1-amino acid dropouts
30626122	CSM-Arg Powder, 10 grams	1-amino acid dropouts
30626123	CSM-Lys Powder, 10 grams	1-amino acid dropouts
30626124	CSM-Arg-His Powder, 10 grams	2-amino acid dropouts
30626125	CSM-Ade-Ura Powder, 10 grams	2-amino acid dropouts
30626126	CSM-Ade-Leu Powder, 10 grams	2-amino acid dropouts
30626127	CSM-Leu-Ura Powder, 10 grams	2-amino acid dropouts
30626128	CSM-Arg-His-Ura Powder, 10 grams	3-amino acid dropouts
30626129	CSM-Arg-His-Lys Powder, 10 grams	3-amino acid dropouts
30626130	CSM-Arg-Leu-Ura Powder, 10 grams	3-amino acid dropouts
30626131	CSM-Arg-Leu Powder, 10 grams	2-amino acid dropouts
30626132	CSM-Leu-Lys Powder, 10 grams	2-amino acid dropouts
30626133	CSM-His-Leu-Lys Powder, 10 grams	3-amino acid dropouts
30626134	CSM-Arg-His-Leu-Trp-Ura Powder, 10 grams	5-amino acid dropouts

For available sizes & current prices, visit: www.plantmedia.com


SKU #	Product Name	Product Group
30626135	CSM-Ade-Arg-His-Leu-Trp Powder, 10 grams	5-amino acid dropouts
30626136	CSM-His-Tyr Powder, 10 grams	2-amino acid dropouts
30626137	CSM-Lys-Trp Powder, 10 grams	2-amino acid dropouts
30626138	CSM-Ade-His-Ura Powder, 10 grams	3-amino acid dropouts
30626139	CSM-Arg-Leu-Trp Powder, 10 grams	3-amino acid dropouts
30626140	CSM-Arg-Trp Powder, 10 grams	2-amino acid dropouts
30626141	CSM-Arg-Trp-Ura Powder, 10 grams	3-amino acid dropouts
30626142	CSM-Leu-Lys-Trp Powder, 10 grams	3-amino acid dropouts
30626143	CSM-Ade-Leu-Trp-Ura Powder, 10 grams	4-amino acid dropouts
30626144	CSM-Ade-His-Leu-Trp-Ura Powder, 10 grams	5-amino acid dropouts
30626145	CSM-His-Leu-Lys-Trp-Ura Powder, 10 grams	5-amino acid dropouts
30626146	CSM-Ade-His-Leu-Ura Powder, 10 grams	4-amino acid dropouts
30626147	CSM-His-Lys-Trp-Ura Powder, 10 grams	4-amino acid dropouts
30626148	CSM-Ade-Trp-Ura Powder, 10 grams	3-amino acid dropouts
30626149	CSM-His-Lys Powder, 10 grams	2-amino acid dropouts
30626150	CSM-His-Leu-Lys-Trp Powder, 10 grams	4-amino acid dropouts
30626151	CSM-Ade-Lys Powder, 10 grams	2-amino acid dropouts
30626152	CSM-Ade-Leu-Ura Powder, 10 grams	3-amino acid dropouts
30626153	CSM-Ade-His-Leu-Lys Powder, 10 grams	4-amino acid dropouts
30626154	CSM-Ade-His-Leu-Lys-Trp-Ura Powder, 10 grams	6-amino acid dropouts
30626155	CSM-Ade-His-Leu-Lys-Trp Powder, 10 grams	5-amino acid dropouts
30626156	CSM-Ade-His Powder, 10 grams	2-amino acid dropouts
30626157	CSM-Ade-His-Trp Powder, 10 grams	3-amino acid dropouts
30626158	CSM-Ade-Arg-His-Ura Powder, 10 grams	4-amino acid dropouts
30626159	CSM-Arg-Leu-Lys-Ura Powder, 10 grams	4-amino acid dropouts
30626160	CSM-Arg-His-Lys-Ura Powder, 10 grams	4-amino acid dropouts
30626161	CSM-His-Lys-Ura Powder, 10 grams	3-amino acid dropouts
30626162	CSM-His-Lys-Trp Powder, 10 grams	3-amino acid dropouts
30626163	CSM-Lys-Trp-Ura Powder, 10 grams	3-amino acid dropouts
30626164	CSM-Leu-Lys-Trp-Ura Powder, 10 grams	4-amino acid dropouts
30626165	CSM-Ade-Arg-His-Leu-Trp-Ura, 10 grams	6-amino acid dropouts
30626166	CSM-Leu-Thr-Trp-Ura Powder, 10 grams	4-amino acid dropouts
30626167	CSM-Leu-Thr-Trp Powder, 10 grams	3-amino acid dropouts
30626168	CSM-Arg-Lys Powder, 10 grams	2-amino acid dropouts
30626169	CSM-Arg-Leu-Lys Powder, 10 grams	3-amino acid dropouts
30626170	CSM-Arg-His-Leu-Lys Powder, 10 grams	4-amino acid dropouts
30626171	CSM-Arg-His-Leu-Lys-Ura Powder, 10 grams	5-amino acid dropouts
30626172	YNB Powder, 250 grams	Yeast Nitrogen Base
30626173	YNB+Nitrogen Powder, 500 grams	Yeast Nitrogen Base
30626174	DOB Powder, 500 grams	Dropout Bases
30626175	DOB 2% Galactose Powder, 250 grams	Dropout Bases
30626176	DOB 2% Raffinose Powder, 250 grams	Dropout Bases

Call 1.888.BIOPUS to order.

SKU #	Product Name	Product Group
30626177	DOB 2% Galactose 1% Raffinose Powder, 250 grams	Dropout Bases
30626178	SD Powder, 500 grams	SD Complete Formulations
30626179	SDA Powder, 500 grams	SD Complete Formulations
30626180	SD-His-Leu-Trp Powder, 500 grams	SD Complete Formulations
30626181	SD-His-Leu-Ura Powder, 500 grams	SD Complete Formulations
30626182	SD-Leu-Trp-Ura Powder, 500 grams	SD Complete Formulations
30626183	SD-His-Leu-Trp-Ura Powder, 500 grams	SD Complete Formulations
30626184	SD-His-Trp-Ura Powder, 500 grams	SD Complete Formulations
30626185	SD-Ade-His-Leu-Trp Powder, 500 grams	SD Complete Formulations
30626186	YPD Powder, 1 kg	YPD
30626187	YPD Agar Slant, 2.5 ml, Box of 25	YPD
30626188	YPAD Powder, 1 kg	YPAD
30626189	EMM Powder, 500 grams	EMM
30626190	YES 50 Powder, 500 grams	YE/YES
30626191	YES 225 Powder, 500 grams	YE/YES
30626192	EMM Low Glucose Powder, 500 grams	EMM
30626193	EMM-Glucose Powder, 500 grams	EMM
30626194	EMM-Nitrogen Powder, 500 grams	EMM
30626195	SP Complete Supplements, 10 grams	0-amino acid dropouts
30626196	SP-Leu-Ura Supplements, 10 grams	2-amino acid dropouts
30626197	SP-Lys Supplements, 10 grams	1-amino acid dropouts
30626198	SP-Leu Supplements Powder, 10 grams	1-amino acid dropouts
30626199	SP-Ade-Leu-Ura Supplements Powder, 10 grams	3-amino acid dropouts



For available sizes & current prices, visit: www.plantmedia.com

SKU #	Product Name	CAS #
Plant Growth Regulators		
		
30631017	(+)- cis, trans-Abscisic acid (ABA)	14375-45-2
42000008	2,4,5-Trichlorophenoxyacetic acid, potassium salt	93-76-5
30631001	2,4-Dichlorophenoxyacetic acid (2,4-D)	94-75-7
30633001	2,4-Dichlorophenoxyacetic acid (2,4-D), Sodium Salt Monohydrate	7084-86-8
30631003	24-Epibrassinolide	78821-43-9
30631004	2-Naphthoxyacetic acid	120-23-0
30631005	6-Benzylaminopurine (6-BAP)	1214-39-7
30631006	6-Benzylaminopurine riboside	4294-16-0
40450004	6-Dimethylallylamino purine riboside	7724-76-7
40450000	6-γ-[Dimethylallylamo]purine (2iP)	2365-40-4
40150000	Adenine (6-Aminopurine)	73-24-5
40100036	Adenine hemisulfate, dihydrate	6509-19-9
40100040	Adenine sulfate	321-30-2
30631016	Alpha-naphthalene acetic acid (NAA)	86-87-3
40150001	Auxindole™	
21520003	Bestatin	58970-76-6
40300002	Coconut Water	
40450016	DICAMBA (3,6-Dichloro-2-methoxybenzoic acid)	1918-00-9
40600056	Fluridone	59756-60-4
30630025	Flurprimidol	56425-91-3
30631009	Gibberellic acid (HPLC)	77-06-5
30631008	Gibberellic acid A3 (GA3)	77-06-5
30631010	Indole-3-acetic acid (IAA)	87-51-4
30631011	Indole-3-butyric acid (IBA)	133-32-4
30631000	Indole-3-Carboxylic acid (ICA)	771-50-6
30631012	Indole-3-propionic acid (IPA)	830-96-6
40950000	Iodoacetamide	144-48-9

Call 1.888.BIOPRO to order.

SKU #	Product Name	CAS #
40950004	Iodoacetic acid, sodium salt	305-53-3
30631013	Jasmonic acid	77026-92-7
30631014	Kinetin (6-Furfuryl-aminopurine)	525-79-1
30633002	Kinetin Solution	525-79-1
30631015	Methyl jasmonate	39924-52-2
41350000	Myo-inositol	87-89-8
41480000	Naphthyl phosphate, sodium salt	81012-89-7
41650004	Paclobutrazol	76738-62-0
30631007	p-Chlorophenoxyacetic acid (CPA)	122-88-3
41650000	Phenylacetic acid (PAA)	103-82-2
41480004	P-Nitrophenyl phosphate	4264-83-9
30631018	Thidiazuron (TDZ)	51707-55-2
30631019	Zeatin	1637-39-4
30631020	Zeatin riboside	6025-53-2
30633003	4-Iodophenoxyacetic Acid	1878-94-0
30633004	1-Triacontanol	593-50-0
30633005	Ethepron	16672-87-0

Plant Media



40100056	Agar, Daishin	9002-18-0
40100060	Agar, Granulated	9002-18-0
40100064	Agar, High Gel Strength	9002-18-0
40100048	Agar, Low Melt	9012-36-3
40100068	Agar, Nobel	9002-18-0
40100072	Agar, Phytoagar	9002-18-0
40100076	Agar, Plant Tissue Culture	9002-18-0
30634002	Agar-Gelzan Blend	

For available sizes & current prices, visit: www.plantmedia.com

SKU #	Product Name	CAS #
30630001	Anderson's Rhododendron Medium	
30630000	Anderson's Rhododendron Medium, w/ Vitamins	
30630002	Bacteria Screening Medium 523	
30630134	Banana AGS Basal Medium	
30630003	Banana Powder	
40200094	Beef Peptone Powder	91079-38-8
30630141	Blaydes Basal Medium	
30630181	BM-1 Terrestrial Orchid Medium	
30630180	BM-1 Terrestrial Orchid Medium, w/ Agar	
30630183	BM-2 Terrestrial Orchid Medium	
30630182	BM-2 Terrestrial Orchid Medium, w/ Agar	
30630156	Cape Sundew/ Venus Flytrap Multiplication Basal Medium	
30630157	Cape Sundew/ Venus Flytrap Pretransplant Basal Medium	
30630089	Carnivorous Plant (Nepenthes) Basal Salt Mixture	
30630158	Carrot Callus Intiation Medium From Pith Tissue	
30630159	Carrot Shoot Development Basal Medium	
30630005	Chee and Pool (C2D) Vitis Medium	
30630006	Chee and Pool (C2D) Vitis Medium, w/ Vitamins	
30630007	Chu (N6) Medium	
30630008	Chu (N6) Medium, w/ Vitamins	
30630009	Chu Vitamin Concentrate 10mL Ampules	
30630010	Chu Vitamin Mix	
30630011	CLC/Ipomoea Basal Medium CP, w/ Vitamins	
30630012	CLC/Ipomoea Basal Medium EP, w/ Vitamins	
40300150	Corn Flour	
30630178	Cornmeal Agar	
30620111	Czapek Dox Agar	
30630014	Czapek Dox Broth	
30630088	DCR Basal Salt Mixture	
30630015	De Greef & Jacobs Medium	
30630016	De Greef & Jacobs Medium, w/ Vitamins	
30630018	De Greef & Jacobs Vitamin Mix	
30630017	De Greef & Jacobs Vitamin Solution (Ampule)	
30630160	DKW Basal Medium, w/ 10 g/L Sucrose	
30630161	DKW Basal Medium, w/ 30 g/L Sucrose	
30630092	DKW Basal Salt Mixture	
30630020	DKW/Juglans Medium, w/ Vitamins	
30630021	DKW/Juglans Vitamin Concentrate (Ampules)	
30630022	DKW/Juglans Vitamin Mix (Concentrate)	

Call 1.888.BIOPUS to order.

SKU #	Product Name	CAS #
30630162	Economou & Read Basal Medium	
30630023	Eriksson (ER) Medium	
30630024	Eriksson (ER) Medium, w/ Vitamins	
30630164	EXS III™ Basal Medium	
30630163	EXS III™ Basal Salt Medium, w/ Adenine	
30630184	Fast Terrestrial Orchid Medium	
30630028	Gamborg B5 Medium	
30630029	Gamborg B5 Medium, w/ Vitamins	
30630026	Gamborg B5 Vitamin Concentration (Ampule)	
30630027	Gamborg B5 Vitamin Mix	
30630090	Gellan Gum Powder	71010-52-1
30634003	Gelzan™ (Gelrite)	71010-52-1
30634004	Gelzan™, Ready-to-use packets	
30630034	Gresshoff & Doy (DBM2) Medium	
30630035	Gresshoff & Doy (DBM2) Medium, w/ Vitamins	
30630033	Gresshoff & Doy (DBM2) Vitamin Concentration (Ampule)	
30630032	Gresshoff & Doy (DBM2) Vitamin Mix	
30630036	Heller Medium	
30630037	Hoagland Complete Medium	
30630038	Hoagland Medium, Nitrogen-Free	
30630165	Hosta Initiation/Multiplication Medium Stage I/II	
30630166	Hosta Multiplication Medium Stage II	
30630167	Hosta Rooting Medium Stage III	
30630196	Ichihashi New Phalaenopsis (NP) Medium	
30630039	Jensen's Medium	
30630040	Kao & Michayluk Medium	
30630041	Knudson C Orchid Medium, Morel Modification	
30630042	Leifert & Waites Sterility Test Medium	
30630043	Lindemann Orchid Medium	
30630171	Linsmaier & Skoog Basal Medium, pH Adjusted and Buffered	
30630168	Linsmaier & Skoog Basal Medium, w/ Sucrose & Agar	
30630169	Linsmaier & Skoog Basal Medium, w/ Sucrose & Agar, pH Adjusted and Buffered	
30630170	Linsmaier & Skoog Basal Medium, w/ Sucrose, pH Adjusted and Buffered	
30630045	Linsmaier & Skoog Medium	
30630046	Linsmaier & Skoog Vitamin Concentration (Ampule)	
30630047	Linsmaier & Skoog Vitamin Mix	
30630048	Litvay Medium	
30630049	Litvay Medium, w/ Vitamins	

For available sizes & current prices, visit: www.plantmedia.com

SKU #	Product Name	CAS #
30630050	Litvay Vitamin Concentrate (Ampule)	
30630051	Litvay Vitamin Mix	
30630185	Malmgren Terrestrial Orchid Medium	
30630186	Malmgren Terrestrial Orchid Medium, w/ Pineapple Powder	
30630052	McCown and Murashige & Skoog (MS) Vitamin Concentration (Ampule)	
30630053	McCown and Murashige & Skoog (MS) Vitamin Mix	
30630056	McCown Woody Plant Basal Salt Mixture	
30630057	McCown Woody Plant Basal Salt Mixture, w/ Vitamins	
30630054	McCown Woody Plant Medium	
30630055	McCown Woody Plant Medium, w/ Vitamins	
30630116	Milieu de Bouturage (S-Medium)	
30630197	Mitra Maintenance/ Replate Medium	
30630083	Murashige & Miller Medium	
30630087	Murashige & Skoog (MS) Basal Salt Concentrate (20X)	
30630070	Murashige & Skoog (MS) Basal Salt Medium, Macronutrients Only	
30630058	Murashige & Skoog (MS) Medium	
30630133	Murashige & Skoog (MS) Medium (1X) Sterile Solution, w/ Vitamins	
30630060	Murashige & Skoog (MS) Medium Plates, w/ Hygromycin	
30630062	Murashige & Skoog (MS) Medium Plates, w/ Kanamycin	
30630069	Murashige & Skoog (MS) Medium Plates, w/ Vitamins	
30630142	Murashige & Skoog (MS) Medium, BC Potato	
30630064	Murashige & Skoog (MS) Medium, Buffered	
30630071	Murashige & Skoog (MS) Medium, Modification No. 1	
30630072	Murashige & Skoog (MS) Medium, Modification No. 1, w/ Vitamins	
30630073	Murashige & Skoog (MS) Medium, Modification No. 2	
30630074	Murashige & Skoog (MS) Medium, Modification No. 2, w/ Vitamins	
30630075	Murashige & Skoog (MS) Medium, Modification No. 3	
30630076	Murashige & Skoog (MS) Medium, Modification No. 3, w/ Vitamins	
30630077	Murashige & Skoog (MS) Medium, Modification No. 4	
30630078	Murashige & Skoog (MS) Medium, Modification No. 5	
30630080	Murashige & Skoog (MS) Medium, w/ Finer & Nagasawa Modification	
30630059	Murashige & Skoog (MS) Medium, w/ Gamborg B5 Vitamins	
30630065	Murashige & Skoog (MS) Medium, w/ Modified Vitamins	
30630066	Murashige & Skoog (MS) Medium, w/ Nitsch Vitamins	
30630136	Murashige & Skoog (MS) Medium, w/ Sucrose and Agar	
30630081	Murashige & Skoog (MS) Medium, w/ Van der Salm Modification	
30630082	Murashige & Skoog (MS) Medium, w/ Van der Salm Modification, w/ Vitamins	
30630067	Murashige & Skoog (MS) Medium, w/ Vitamins	

Call 1.888.BIOPRO to order.

SKU #	Product Name	CAS #
30630068	Murashige & Skoog (MS) Medium, w/ Vitamins, Buffered	
30630143	Murashige & Skoog (MS) Modified Medium, Arabidopsis Culture Medium	
30630144	Murashige & Skoog (MS) Modified Medium, w/ Adenine Hemisulfate	
30630145	Murashige & Skoog (MS) Modified Medium, w/ IAA	
30630146	Murashige & Skoog (MS) Modified Multiplication Medium	
30630086	Murashige & Skoog (MS) Modified Vitamin Powder	
30630079	Murashige & Skoog (MS) Shoot Multiplication B Medium	
30630129	Murashige & Skoog (MS) Stock Solution (10X), Macronutrient Only	
30630128	Murashige & Skoog (MS) Stock Solution (10X), Micronutrient Only	
30630130	Murashige & Skoog (MS) Vitamin Mix Concentrate (Ampule)	
30630147	Murashige African Violet/Gloxinia Multiplication Medium	
30630148	Murashige African Violet/Gloxinia Pretransplant Medium	
30630187	Murashige Cattleya Orchid Multiplication Medium	
30630149	Murashige Fern Multiplication Medium	
30630150	Murashige Gerbera Multiplication Medium	
30630151	Murashige Gerbera Pretansplant Medium	
30630152	Murashige Kalanchoe Multiplication Medium	
30630153	Murashige Kalanchoe Pretransplant Medium	
30630154	Murashige Modified Multiplication Medium, w/ Kinetin	
30630155	MUSA (Banana) Multiplication Medium	
30630172	NB Basal Medium	
30620040	Nematode Growth Medium (NGM)	
30630095	Nitsch Medium	
30630096	Nitsch Medium, w/ Vitamins	
30630097	Nitsch Vitamin Concentrate (Ampule)	
30630098	Nitsch Vitamin Mix	
30630102	NLN Medium	
30630103	NLN Medium, w/ Vitamins	
30630179	Oatmeal Agar	
30630188	Orchid Maintenance Medium, w/ Charcoal and Agar	
30630189	Orchid Maintenance/ Replate Medium, w/ Banana and Charcoal	
30630190	Orchid Maintenance/ Replate Medium, w/ Banana, Charcoal, and Agar	
30630191	Orchid Multiplication Medium	
30630193	Orchid Seed Sowing Medium, w/ Agar	
30630192	Orchid Seed Sowing Medium, w/ Agar & Charcoal	
30630105	Orchimax Medium	
30630104	Orchimax Medium, w/ Charcoal	
30630139	Parker Thompson Fern Basal Salt Mixture	
30630107	Quoirin & Lepoivre Medium	

For available sizes & current prices, visit: www.plantmedia.com

SKU #	Product Name	CAS #
30630108	Quoirin & Lepoivre Medium, w/ Vitamins	
30630109	RDM (Rhizobium defined medium)	
30630173	Rose Initiation Medium	
30630174	Rose Modified Multiplication Basal Medium	
30630110	Rugini Olive Medium	
30630111	Schenk & Hildebrandt Medium	
30630135	Schenk & Hildebrandt Medium, w/ Vitamins	
30630176	Schenk & Hildebrandt Modified Medium, w/ Sucrose	
30630093	Schenk & Hildebrandt Modified Salt Mixture, w/o CaCl3	
30630113	Schenk & Hildebrandt Vitamin Concentration (Amplule)	
30630112	Schenk & Hildebrandt Vitamin Mix	
30630114	Sequestrene™ 138, 6% Iron chelate	
30630115	Sequestrene™ 330, 10% Iron chelate	
30630194	Terrestrial (Cypripedium) Orchid Medium	
30630138	TM4G Basal Salt Mixture	
30630177	TM4G Medium	
30630117	Tomato Powder	
30630118	Vacin & Went Medium	
30630140	Vacin & Went Modified Orchid Basal Salt Mixture	
30630195	Vacin & Went Orchid Medium, w/ Agar and Sucrose	
30630119	Westvaco WV3 Medium	
30630120	Westvaco WV5 Medium	
30630121	Wheat Callus Induction Medium (CM4C)	
30630122	Wheat Regeneration Medium (MMS2C)	
30630123	Wheat Rooting Medium (MMSOC)	
30630124	White's Medium	
30630125	YEB Agrobacterium Growth Medium	
30630126	YPD Agar	
30630127	YT Broth	

Plant Media Chemicals



40100297	Acetosyringone (3',5'-Dimethoxy-4'-Hydroxyacetophenone)	2478-38-8
40100028	Activated charcoal acid washed	

Call 1.888.BIOPRO to order.

SKU #	Product Name	CAS #
40100204	Ammonium Nitrate	6484-52-2
40100208	Ammonium phosphate, dibasic	7783-28-0
40100212	Ammonium phosphate, monobasic	7722-76-1
40110032	Ampicillin, sodium salt	69-52-3
40220048	Borate Buffer, 0.05M, (pH 8.5)	
40300012	Calcium Chloride, anhydrous	10043-52-4
40300151	Calcium Hydrogen Phosphate	7757-93-9
40300032	Calcium nitrate	
40300153	Calcium Phosphate, Tribasic	12167-74-7
40310008	Cefotaxime Sodium	64485-93-4
30900014	Citric Acid (Plant Tissue Culture Grade)	77-92-9
40300124	Cyanocobalamin	68-19-9
40600008	D-Fructose	57-48-7
40700008	D-Glucose, anhydrous	50-99-7
40700012	D-Glucose, monohydrate	5996-10-1
41300001	D-Maltose, monohydrate	6363-53-7
41300012	D-Mannitol	69-65-8
41300016	D-Mannose	3458-28-4
41800000	D-Ribose	50-69-1
40500032	Ergo-calciferol	50-14-6
41900216	Sodium Phosphate, Monobasic, Anhydrous	7558-8--7

Herbicides



30632000	Atrazine	1912-24-9
30632006	Basta (1% solution)	
30632007	Basta (Phosphinothrin/Glufosinate Ammonium)	77182-82-2
40200093	Bromoxynil	1689-84-5
30632001	Chlorsulfuron	64902-72-3
30632002	DL-Phosphinothrin, PPT	77182-82-2

For available sizes & current prices, visit: www.plantmedia.com

SKU #	Product Name	CAS #
30630001	Anderson's Rhododendron Medium	
30630000	Anderson's Rhododendron Medium, w/ Vitamins	
30630002	Bacteria Screening Medium 523	
30630134	Banana AGS Basal Medium	
30630003	Banana Powder	
40200094	Beef Peptone Powder	91079-38-8
30630141	Blaydes Basal Medium	
30630181	BM-1 Terrestrial Orchid Medium	
30630180	BM-1 Terrestrial Orchid Medium, w/ Agar	
30630183	BM-2 Terrestrial Orchid Medium	
30630182	BM-2 Terrestrial Orchid Medium, w/ Agar	
30630156	Cape Sundew/ Venus Flytrap Multiplication Basal Medium	
30630157	Cape Sundew/ Venus Flytrap Pretransplant Basal Medium	
30630089	Carnivorous Plant (Nepenthes) Basal Salt Mixture	
30630158	Carrot Callus Initiation Medium From Pith Tissue	
30630159	Carrot Shoot Development Basal Medium	
30630005	Chee and Pool (C2D) Vitis Medium	
30630006	Chee and Pool (C2D) Vitis Medium, w/ Vitamins	
30630007	Chu (N6) Medium	
30630008	Chu (N6) Medium, w/ Vitamins	
30630009	Chu Vitamin Concentrate 10mL Ampules	
30630010	Chu Vitamin Mix	
30630011	CLC/Ipomoea Basal Medium CP, w/ Vitamins	
30630012	CLC/Ipomoea Basal Medium EP, w/ Vitamins	
40300150	Corn Flour	
30630178	Cornmeal Agar	
30620111	Czapek Dox Agar	
30630014	Czapek Dox Broth	
30630088	DCR Basal Salt Mixture	
30630015	De Greef & Jacobs Medium	
30630016	De Greef & Jacobs Medium, w/ Vitamins	
30630018	De Greef & Jacobs Vitamin Mix	
30630017	De Greef & Jacobs Vitamin Solution (Ampule)	
30630160	DKW Basal Medium, w/ 10 g/L Sucrose	
30630161	DKW Basal Medium, w/ 30 g/L Sucrose	
30630092	DKW Basal Salt Mixture	
30630020	DKW/Juglans Medium, w/ Vitamins	
30630021	DKW/Juglans Vitamin Concentrate (Ampules)	
30630022	DKW/Juglans Vitamin Mix (Concentrate)	

Call 1.888.BIOPUS to order.

SKU #	Product Name	CAS #
40710012	Gamma-cyclodextrin, hydrate	91464-90-3
30632003	Glyphosate	1071-83-6
41300010	Metribuzin	21087-64-9
41650008	Picloram	1918-02-1
Plant Antibiotics		
40110028	Ampicillin HP	
42010012	Timentin (Ticarcillin/Clavunate, 15:1)	4697-14-76/ 61177-45-5
40310000	Carbenicillin disodium salt	4800-94-6
41110000	Kanamycin sulfate	25389-94-0
40710026	G418 Sulfate	108321-42-2
30630004	BASTA	

***Custom Media
Formulations Available***



Call 1.888.BIOPPLUS to order.

Comparison of Different MS (Murashige & Skoog) Based Media

SKU#	Product Name	Vitamins	Hormones					Sucrose g/L	Agar	Other Ingredients
			IAA mg/L	NAA mg/L	6-BA mg/L	2iP mg/L	Kinetin mg/L			
30630058	Murashige & Skoog (MS) Medium									Macro & Micronutrients only
30630067	Murashige & Skoog (MS) Medium w/ Vitamins	MS								
30630136	Murashige & Skoog (MS) Medium w/ Sucrose & Agar	MS						30	7	
30630086	Murashige & Skoog (MS) Modified Vitamin Mix	Modified MS								Thiamine (1000mg)
30630059	Murashige & Skoog (MS) Medium with Gamborg's Vitamins	GB5								
30630066	Murashige & Skoog (MS) Medium w/ Nitsch Vitamins	Nitsch								
30630064	Murashige & Skoog (MS) Medium, Buffered	MS								Contains MES Buffer
30630070	Murashige & Skoog (MS) Medium, Macro salt only									Macro salts only
30630202	Murashige & Skoog (MS) Micro									Micro salts only
30630133	Murashige & Skoog (MS) Medium, 1X Sterile solution	Ms						30		Macro & micro salts, vitamins, sucrose
30630129	Murashige & Skoog (MS) 10X Stock solution, Macro only									Macronutrients only (liquid)
30630128	Murashige & Skoog (MS) Micro									Micronutrient salts only (liquid)
30630198	Murashige & Skoog (MS) Vitamin Powder Mix	MS								
30630130	Murashige & Skoog (MS) Vitamin Mix	MS								1000X, liquid
30630071	Murashige & Skoog (MS) Medium, Modification 1						0.1			1/2 Micronutrients
30630072	Murashige & Skoog (MS) Medium, Modification 1	MS								1/2 Micronutrients
30630073	Murashige & Skoog (MS) Medium, Modification 2									3/4 Macronutrients
30630074	Murashige & Skoog (MS) Medium, Modification 2	MS								3/4 Macronutrients
30630075	Murashige & Skoog (MS) Medium, Modification 3									1/2 Ammonium Nitrate and Potassium Nitrate
30630076	Murashige & Skoog (MS) Medium, Modification 3	MS								1/2 Ammonium Nitrate and Potassium Nitrate
30630077	Murashige & Skoog (MS) Medium, Modification 4									Ammonium nitrate free
30630078	Murashige & Skoog (MS) Medium, Modification 5									Ammonium nitrate replaced by Sodium nitrate
30630144	Murashige & Skoog (MS) Modified medium	LS*								Adenine Hemisulfate (80mg/L)
30630201	Murashige & Skoog (MS) Modified medium w/ o Potassium Phosphate									w/o Potassium Phosphate

Gamborg (GB5) vitamins: Myo-Inositol (100mg/L), Nicotinic acid (1.0mg/L), Pyridoxine HCl (1mg/L) and Thiamine HCL (10mg/L)

LS vitamins: Myo-Inositol (100mg/L) & Thamine HCL (0.40 mg/L)

MS vitamins: Glycine (2mg/L), Myo-Inositol (100mg/L), Nicotinic acid (0.5mg/L), Pyridoxine HCL (0.5mg/L) and Thamine HCL (0.1mg/L)

Nitsch vitamins: Biotin (0.05mg/L), Folic acid (0.5mg/L), Glycine (2mg/L), Myo-Inositol (100mg/L), Nicotinic acid (5mg/L), Pyridoxine HCL (0.5mg/L) & Thiamine HCL (0.5mg/L)

* Linsmaier & Skoog (LS) medium contains the same macro and micronutrients as Murashige & Skoog medium differing only in its vitamin content.

For available sizes & current prices, visit: www.plantmedia.com

Comparison of Different MS (Murashige & Skoog) Based Media

SKU#	Product Name	Vitamins	Hormones					Sucrose g/L	Agar	Other Ingredients
			IAA mg/L	NAA mg/L	6-BA mg/L	2iP mg/L	Kinetin mg/L			
30630145	Murashige & Skoog (MS) Modified Medium w/ IAA	LS*	0.3					30		Sodium Phosphate (148mg/L) Adenine Hemisulfate (80mg/L)
30630199	Murashige & Skoog (MS) Medium, Modified with FeNaEDTA	MS								FeNaEDTA replaces FeSO4 and NaEDTA
30630200	Murashige & Skoog (MS) Medium w/o Nitrogen									W/o ammonium nitrate and Potassium nitrate
30630203	Murashige & Skoog (MS) Modified Medium									1/2 macro and Micro
30630204	Murashige & Skoog (MS) Modified medium									1/2 NH4NO3, KNO3, CaCl2
30630207	Murashige & Skoog (MS) Modified medium with NAA and BA			0.1	0.1			30		
30630206	Murashige & Skoog / Gamborg (MG) Medium									Modified Murashige/Gamborg salts
30630205	Murashige & Skoog (MS) Modified medium									w/o Nitrogen, Potassium and Phosphates
30630079	Murashige & Skoog (MS) Shoot Multiplication Medium B	LS*								Sodium Phosphate (170mg/L), Adenine Hemisulfate (80mg/L) & 0.4mg/L of Thiamine
30630146	Murashige & Skoog (MS) Modified Multiplication Medium C	LS*		0.1			1	30		
30630154	Murashige & Skoog (MS) Modified Multiplication Medium	LS*	0.3				1			
30630149	Murashige Modified Fern Multiplication Medium	LS*		0.1			2			Sodium Phosphate (255mg/L)
30630150	Murashige Modified Gerbera Multiplication Medium	MS	0.5				10			Sodium Phosphate (85mg/L), Adenine sulfate (80mg/L) & L-Tyrosine (100mg/L)
30630151	Murashige Modified Gerbera Pretransplant Medium	MS	10							Sodium Phosphate (85mg/L) & L-Tyrosine (100mg/L)
30630152	Murashige Modified Kalanchoe Multiplication Medium	LS*				3				
306300080	Murashige & Skoog (MS) Medium, Finer & Nagasawa Modification	MS								1.6X KNO3, 0.5x NH4NO3 and FeNaEDTA in place of FeSO4 and NaEDTA
30630142	Murashige & Skoog (MS) BC Potato Medium	Modified MS					0.04			
30630147	Murashige Modified African Violet/ Gloxinia Multiplication Medium	LS*	2				2			Sodium Phosphate (170mg/L), Adenine Hemisulfate (80mg/L)
30630148	Murashige Modified African Violet/ Gloxinia Pretransplant Medium	LS*	1							
30630187	Murashige & Skoog (MS) Cattleya Orchid Multiplication Medium	MS	0.3	1.75	1.75			20		
30630081	Murashige & Skoog (MS) Medium, Van Der Salm Modification									FeNaEDDHA replaces FeNaEDTA
30630082	Murashige & Skoog (MS) Cattleya Orchid Multiplication Medium	MS								FeNaEDDHA replaces FeNaEDTA
30630045	Linsmair & Skoog medium	LS*								
30630168	Linsmair & Skoog Basal Medium w/ Sucrose & Agar	LS*						30	7	

All media are used with 0.7% phytoagar. Catalog # 40100072

Call 1.888.BIOPRO to order.



plantmedia
a division of bioWORLD



1.888.bio.PLUS
www.bio-world.com